

Optimizing Protocols for Automated Multistep Purification of (histidine)₆-tagged Protein using ÄKTExpress

A. Sjöberg¹, N. Michot², E. Linde¹, M. Nilsson¹, R. Bhikhabhai¹, I. Salomonsson¹, L. Hedkvist¹, M. Galin¹, P. Liljedahl¹, C. Markeland-Johansson¹, P. Eklund¹, T. Frigård¹, N. Pettersson¹, J. Dumas² and J.A. Simon¹.

¹ GE Healthcare, SE-751 84 Uppsala, Sweden

² Aventis Pharma, Paris Research Center, Lead Discovery and Technologies/Protein Production, 94 403 VITRY sur SEINE, France

Introduction

ÄKTExpress™ is a chromatography system for automated multistep purification of large numbers of (histidine)₆- and Glutathione S-transferase (GST)-tagged proteins. Pre-optimized methods are easily created using the method wizard in UNICORN™ control software and can be further optimized using a special feature called Advanced Zone.

To increase the final purity of a target protein, parameters for each individual step can be changed. In this paper, two different optimization approaches purifying a (histidine)₆-tagged kinase are described.

- Increased imidazole concentration in the wash step prior to affinity elution
- Shallow gradient during ion exchange elution

ÄKTExpress

- Automated multistep purification of affinity-tagged proteins
- Method wizard for easy creation of purification protocols
- Intelligent peak detection and collection in intermediate steps
- Optional on-column tag cleavage
- Up to four samples can be purified per module
- Up to twelve modules in parallel can be controlled from one computer



Fig 1. ÄKTExpress system containing 2 modules.

Conclusions

The Advanced Zone allows further optimization of ÄKTExpress purification protocols.

For a three-step protocol (Affinity – Desalting – Ion exchange),

- purifications performed with default values yielded pure test protein
- a wash with raised imidazole concentration, prior to elution in the affinity step, resulted in increased purity of the test protein
- higher resolution was obtained using a shallower gradient in the ion exchange step



All protocols start with affinity chromatography followed by different combinations of desalting, ion exchange chromatography, and gel filtration. The largest peak from each step is transferred to the next column.

Fig 2. Schematic chromatogram of an automated four-step protein purification protocol.




Fig 2. Schematic chromatogram of an automated four-step protein purification protocol.

Table 1. Descriptions of available multistep protocols supported by ÄKTApres

Available multistep protocols	Effects from additional chromatographic steps
AC-DS	Buffer exchange
AC-GF	Separation from undesired aggregates and contaminants
AC-DS-IEX	Separation from other isoforms (e.g. heterogenously phosphorylated or glycosylated proteins)
AC-DS-IEX-DS	Separation from other isoforms on IEX and buffer exchange on DS
AC-DS-IEX-GF	Separation from other isoforms on IEX and removal of undesired aggregates and contaminants on GF

Material and methods

Sample

A (histidine)₆-tagged kinase (M_r 42.4x10³, pI 5.75) expressed in *E. coli* was used in this study.

Buffers

AC binding buffer:	50 mM Tris-HCl, 500 mM NaCl, 5 mM imidazole, 1 mM DTT, 10% glycerol, pH 8
AC elution buffer:	50 mM Tris-HCl, 500 mM NaCl, 500 mM imidazole, 1 mM DTT, 10% glycerol, pH 8
DS and IEX binding buffer:	50 mM Tris-HCl, 25 mM NaCl, 1 mM DTT, 10% glycerol, pH 7.5
IEX elution buffer:	50 mM Tris-HCl, 1 M NaCl, 1 mM DTT, 10% glycerol, pH 7.5
GF buffer:	50 mM Tris-HCl, 150 mM NaCl, 1 mM DTT, 10% glycerol, pH 7.5

Columns

AC:	HisTrap™ HP, 1 ml
DS:	HiPrep™ 26/10 Desalting
IEX:	Mono Q™ 5/50 GL, 1 ml
GF:	HiLoad™ 16/60 Superdex™ 75 pg

Analysis

Purity was analyzed on Coomassie™ stained ExcelGel™ 8–18 % SDS polyacrylamide gels under reduced conditions. Approximately 7.5 µg of protein was loaded per lane.

Advanced Zone

When creating ÄKTExpress methods, it is possible to view or edit parameters from the Advanced Zone in the method wizard. A selection of the parts of the method to be changed are made on the first page in Advanced Zone (Fig 3). In the following pages, the variables are shown and can be edited in a straightforward way.

There are many variables that can be changed in Advanced Zone; some of these are described below.

Additional affinity wash: By using harsher wash conditions (e.g. a raised concentration of imidazole) weak and nonspecifically bound proteins are removed. This can increase final purity.

Ion exchange gradient: To increase the resolution of an ion exchange column, elution flow rate can be decreased or the gradient parameters (length and target concentration) can be altered.

Peak detection: For peak detection, values are column dependent and can be altered when default settings do not give desired results.





Fig 3. Two pages from the method wizard's Advanced Zone.

Table 2. Examples of variables used in a four-step protocol

Technique	Variable	Default setting
Affinity (AC)	Additional wash Parameters for peak detection	0% B during 20 CV Column dependent
Desalting (DS)	Parameters for peak detection	Column dependent
Ion exchange (IEX)	Flow rates Gradient settings Parameters for peak detection	Column dependent 0–50% B during 20 CV Column dependent
Gel filtration (GF)	Parameters for peak fractionation	Column dependent

Optimization of the affinity and ion exchange steps


AC-DS-IEX default



Lane 1. Low Molecular Weight marker
 Lane 2. Start sample
 Lane 3. Flow through
 Lane 4-9. See corresponding fractions in chromatogram
 Lane 10. Low Molecular Weight marker

Fig 4.


AC-DS-IEX with additional wash



Lane 1. Low Molecular Weight marker
 Lane 2-9. See corresponding fractions in chromatogram

Fig 5.

AC-DS-IEX with 80 CV IEX gradient



Lane 1-9. See corresponding fractions in chromatogram
 Lane 10. Low Molecular Weight marker

Fig 6.

Discussion

- An additional wash prior affinity elution removed the weakly bound M_r 66x10³ protein and low molecular weight contaminants (compare the first peak in the chromatograms of Fig 4 and 5, as well as SDS-PAGE lane 4-6 in Fig 4 with lane 2-4 in Fig 5).
- The resolution was slightly increased using a shallower ion exchange elution gradient compared to when using the default gradient setting (Fig 6).

- By adding a gel filtration column as the fourth step, the protein was conditioned into final buffer and size homogeneity guaranteed. In this study however, the purity achieved after both three- and four-step purifications was very similar (results not shown).
- Another way to increase final purity would be to overload the affinity column. By this approach, strong binders such as Histidine-tagged proteins, will out-compete weak binding compounds during sample loading (results not shown).