

List of chemical resistance of materials in contact with liquid

ANALYSETTE 22 NeXT & ANALYSETTE 28 ImageSizer

LASER PARTICLE SIZER & PARTICLE SIZER

This chart information is based on FRITSCH's knowledge and experience and are intended as a guide and not as a guarantee. We cannot assure the perfect performance of the materials as that depends on many factors as working pressure, pressure picks, fluid, ambient temperature, fluid temperature, concentration, duration of exposure, etc...



Fritsch GmbH
Milling and Sizing
Industriestrasse 8
D - 55743 Idar-Oberstein
Telephone: +49 (0) 6784/ 70-0
Email: info@fritsch.de
Internet: www.fritsch.de

Table of contents

1	Materials in contact with liquid.....	4
2	Abbreviations.....	5
3	List of chemical resistance.....	7
3.1	KEY.....	7
3.2	Note.....	7
3.3	Table 1.....	7
3.4	Table 2.....	97
3.5	Table 3.....	105

Materials in contact with liquid

1 Materials in contact with liquid

Wet dispersion units	Materials used in the sample circuit:
Wet dispersion unit (22.9200.00 & 28.2500.00)	stainless steel ¹ , PTFE (Teflon) ² , N-BK7 glass, Viton ^{®2} , Viton [®] Extreme ³ , POM ^a , Brass, Silicone hoses (LEZ-SIL 60 PFT) ^o
"Conversion Kit Extended" for extreme chemical resistance (22.9287.00 & 22.9288.00)	FFKM (Kalrez [®]) ³ : replaces Viton ^{®2} & Viton [®] Extreme ³ ; TPE (LEZ-SAN Santoprene [®]) ³ : replaces Silicone hoses (LEZ-SIL 60 PFT)

¹ See Table 1 and Table 2

² See Table 2

³ See Table 1

^o See Table 2 (Silicone (peroxide/ platinum- cured))

^a See Table 3

2 Abbreviations

ASTM D1418 Abbreviation	Chemical name	Trade names	PPE Prefix
ACM	Polyacrylate (Polyacrylic)	Europrene-AR®, Krynac®, Nipolar®, Hycar®, Cyanacryl®, Noxtite-PA®	P
AEM	Ethylene Acrylic	Vamac® W	W
AQU / EU	Polyester Urethane Polyether Urethane	Vibrathane®, Pellethane®, Urepan®, Elastothane®, Adiprene®, Desmopan®, Vulkollan®, Estane®	Q
CR	Chloroprene	Neoprene®, Skyprene®, Butclor®, Baypren®, Denka®	C
CSM	Chlorosulphonylpolyethylene (Chlorosulphonated polyethylene)	Hypalon®, Noralon®	L
ECO	Epichlorohydrin	Hydrin®, Herclor®, Gechron®	H
EPD / EPDM	Ethylene-Propylene Copolymer, Ethylene-Propylene Terpolymer	Buna-AP®, Dutral®, Epsyn®, Nordel®, Epcar®, Keltan®, Royalene®, Polysar-EDM®	E
FEP / PFA	Fluoroethylene Propylene, Perfluoroalkoxy	Teflon FEP®, Neoflon FEP®, Teflon PFA®, Neoflon PFA®	D
FEPM	See TFE/P	Aflas®, Dyneon BRE®; Viton Extreme®	A
FFKM/FFPM*	Perfluoroelastomer	Perlast®, Kalrez®, Chemraz®	G
FKM/FPM*	Fluoroelastomer	Viton®, Dai-el®, Fluorel®, Tecnoflon®	V
FVMQ	Fluorosilicone	Silastic, LS®, FSE®	F
NBR / XNBR	Nitrile Butadiene, Carboxylated Nitrile	Breon, Butakon®, Chemigum®, Hycar®, Buna-N®, Butacril®, Paracil®, Perbunan®, Krynac®, Europrene-N®, Nipol®	N
HNBR	Hydrogenated Nitrile	Therban®, Tornac®, Zetpol®	Z
IIR	Butyl Rubber	Exxon Butyl®, Polysar Butyl®, Esso Butyl®	J
IR	Isoprene	Shell Isoprene Rubber®	R
NR	Natural Rubber	Natsyn®	R
POM	Polyoxymethylene		

Abbreviations

ASTM D1418 Abbreviation	Chemical name	Trade names	PPE Prefix
SBR	Styrene Butadine (Buna-S)	Cariflex S [®] , Plioflex [®] , Euro-prene [®] , Pliolite [®] , Buna Huls [®] , Carom [®] , Solprene [®]	B
TFE/P	Tetrafluoroethylene/Propylene	Aflas [®] , Flourel II [®] , Fluoraz [®]	A
TPE	Thermoplastisches Elastomer	LEZ-SAN Santoprene [®]	
VMQ / PVMQ	Silicone	Silastic [®] , Siloprene [®] , Rhodorsil [®] , Silplus [®]	S

3 List of chemical resistance

3.1 KEY

A - No effect; little noticeable change

B - Minor effect; slight corrosion or discoloration

C - Moderate effect; not recommended for continuous use; softening, loss of strength, swelling and/or shrinkage

D - Severe effect; not recommended for use; severe softening, swelling and/or shrinkage

3.2 Note

The CHEMICAL PROPERTIES of N-BK7 glass see: https://www.schott.com/d/advanced_optics/c36214d9-13c4-468c-bf40-8d438b89f532/1.30/schott-optical-glass-pocket-catalog-jan-2020-row.pdf

3.3 Table 1

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Abietic Acid		A	A		
Acetaldehyde	D	A	A	D	A
Acetamide	C	C	A	A	A
Acetanilide	C	A			
Acetic Acid, 30%		A	A	A	A
Acetic Acid, 5%	A	A	A	A	A
Acetic Acid, Glacial	B	A	A	A	A
Acetic Acid, Hot, High Pressure	D	B	A		A
Acetic Anhydride	D	C	A	A	A
Acetoacetic Acid	C	C	A		A
Acetone	D	B	A	A	A
Acetone Cyano-hydrin	C	B	A	A	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Acetonitrile	A	B	A	A	C
Acetophenetidine	A		A		
Acetophenone	D	C	A	B	
Acetotoluidide	A		A		
Acetyl Acetone	D	A	A	A	
Acetyl Bromide	A		A	D	
Acetyl Chloride	A	A	A	D	A
Acetylene	A	A	A	A	A
Acetylene Tetra-bromide	A	A	A		
Acetylene Tetra-chloride	A	A	A	D	
Acetylsalicylic Acid	A	A	A		
Acids, Non-organic		A	A		
Acids, Organic		A	A		
Aconitic Acid			A		
Acridine			A		
Acrolein	C	C	A	A	
Acrylic Acid	A	B		A	
Acrylonitrile	C	C	B	D	
Adipic Acid		A	B	B	
Aero Lubriplate	A	A	A		
Aero Shell 17 Grease	A	A	A		
Aero Shell 750	A	A	A		
Aero Shell 7A Grease	A	A	A		
Aero Shell IAC	A	A	A		
Aerosafe 2300	D	C	B		
Aerosafe 2300W	D	C	B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Aerozene 50 (50% Hydrazine 50% UDMH)	D	C	B		
Air, Below 200°F	A	A	A		
Air, 200 - 300°F	A		A	D	
Air, 300 - 400°F	A		A		
Air, 400 - 500°F	C		B		
Aliphatic Dicarboxylic Acid	A		A	D	
Alkanes (Paraffin Hydrocarbons)	A	A	A		
Alkanesulfonic Acid	A		A		
Alkazene	B	A	A		
Alkenes (Olefin Hydrocarbons)	A	A	A		
Alkyl Acetone	C	C	A		
Alkyl Alcohol	A	A	A		
Alkyl Amine	A	B	B		
Alkyl Aryl Sulfonates	A		A		
Alkyl Aryl Sulfonics	A		A		
Alkyl Benzene	A	A	A		
Alkyl Chloride	A	A	A		
Alkyl Sulfide*	A		A		
Alkylnaphthalene Sulfonic Acid	A		A		
Allyl Chloride	A		A	C	
Allylidene Diacetate	C		B		
Alpha Picoline	C		A		
Aluminum Acetate	D	C	A	A	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Aluminum Bromide	A	A	A	A	
Aluminum Chlorate	C	A	A	A	
Aluminum Chloride	A	A	A	A	B
Aluminum Ethylate			A	A	
Aluminum Fluoride	A		A	A	C
Aluminum Fluorosilicate*			A	A	
Aluminum Formate	C	C	A	A	
Aluminum Hydroxide	B	A	A	A	B
Aluminum Linoleate	A		A	A	
Aluminum Nitrate	A	A	A	A	A
Aluminum Oxalate	C		A		
Aluminum Phosphate	A	A	A		
Aluminum Potassium Sulfate	C	A	A		B
Aluminum Salts	A	A	A	A	
Aluminum Sodium Sulfate	C	A	A		
Aluminum Sulfate	A	A	A	A	B
Alums-NH3 -Cr -K	D	A	A		
Ambrex 33 (Mobil)	A	A	A		
Ambrex 830 (Mobil)	A	A	A		
Amines-Mixed	D	C	A	C	A
Aminoanthraquinone			A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Aminoazobenzene			A	A	
Aminobenzene Sulfonic Acid			A		
Aminobenzoic Acid			A		
Aminopyridine			A		
Aminosalicylic Acid			A		
Ammonia (Anhydrous)	D	B	A		A
Ammonia and Lithium Metal in Solution	D				
Ammonia, Gas, Cold	D		A		
Ammonia, Gas, Hot	D		A		
Ammonia, Liquid (Anhydrous)	D				A
Ammonium Acetate	C	C	A	A	
Ammonium Arsenate	C		A		
Ammonium Benzoate	C		A		
Ammonium Bicarbonate	C	C	A		
Ammonium Bisulfite	C		A		
Ammonium Bromide	A		A		
Ammonium Carbamate	C		A		
Ammonium Carbonate	A		A	A	A
Ammonium Chloride, 2N	A	A	A	A	C

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ammonium Citrate	C		A		
Ammonium Dichromate	C		A		
Ammonium Diphosphate	C	A	A		
Ammonium Fluoride	A	A	B		
Ammonium Fluorosilicate*			A		
Ammonium Formate	C		A		
Ammonium Hydroxide, 3 Molar	C	A	B	A	A
Ammonium Hydroxide, Concentrated	D	A	B	A	A
Ammonium Iodide	A		A		
Ammonium Lactate	C		A		
Ammonium Metaphosphate	C		A		
Ammonium Molybdenate*	C		A		
Ammonium Nitrate, 2N		A	A	A	A
Ammonium Nitrite		A	A		A
Ammonium Oxalate	C		A		A
Ammonium Perchlorate	C		A		
Ammonium Perchloride			A		
Ammonium Persulfate 10%		A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ammonium Per-sulfate Solution		A	A		
Ammonium Phosphate	D	A	A		
Ammonium Phosphate, Dibasic		A	A		A
Ammonium Phosphate, Mono-Basic		A	A		A
Ammonium Phosphate, Tri-basic		A	A		A
Ammonium Phosphite	C		A		
Ammonium Pic-rate	C		A		
Ammonium Poly-sulfide	C		A		
Ammonium Sali-cylate	C		A		
Ammonium Salts	C	B	A	A	
Ammonium Sul-famate	C	C	A		
Ammonium Sul-fate	D	A	A	A	B
Ammonium Sul-fate Nitrate	D		A		
Ammonium Sul-fide	D	A	A		
Ammonium Sul-fite	C	A	A	A	
Ammonium Thio-cyanate	C	A	A		
Ammonium Thio-glycolate	C		A		
Ammonium Thio-sulfate	C	A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ammonium Tungstate	C		A		
Ammonium Valerate	C		A		
Amyl Acetate	D	B	A	D	A
Amyl Alcohol	B	A	A	A	A
Amyl Borate	A	A	A		
Amyl Butyrate	A	A	A		
Amyl Chloride	A	A	A	D	B
Amyl Chloronaphthalene	A	A	A		
Amyl Cinnamic Aldehyde	A	C	A		
Amyl Laurate	A	A	A		
Amyl Mercaptan	A	A	A		
Amyl Naphthalene	A	A	A		
Amyl Nitrate	C		A		
Amyl Phenol		A	A		
Amyl Propionate	A	A	A		
Anderol, L- 826 (di-ester)	A	A	A		
Anderol, L- 829 (di-ester)	A	A	A		
Anderol, L-774 (di-ester)	A	A	A		
ANG-25 (Di-ester Base) (TG749)	A	A			
ANG-25 (Glyceral Ester)	A	A			
Aniline	C	A	A	A	A
Aniline Dyes	B	A	A	A	
Aniline Hydrochloride	B	A	A		
Aniline Oil	C	B	A	B	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Aniline Sulfate	C		A		
Aniline Sulfite	C		A		
Animal Fats	A	A	A	B	
Animal Oil (Lard Oil)	A	A	A		
Anisole		A	A		
Anisoyl Chloride			A		
AN-O-3 Grade M	A				
AN-O-366	A				
AN-O-6	A				
Ansul Ether 161 or 181	D	C	A		
Anthracene	A	A	A		
Anthranilic Acid		A	A		
Anthraquinone		A	A		
Anti-freeze Solutions	C	A	A		
Antimony Chloride	A	A	A		
Antimony Pentachloride	A	A	A	A	
Antimony Pentafluoride		A	B		
Antimony Sulfate			A		
Antimony Tribromide	A		A		
Antimony Trichloride	A		A		D
Antimony Trifluoride	A		B		
Antimony Trioxide	A		A		
AN-VV-O-366b Hydr. Fluid	A	A	A		
Aqua Regia	B	A	A	D	D

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Arachidic Acid		A	A		
Argon	A	A	A	A	
Aroclor, 1248	A	A	A		
Aroclor, 1254	A	A	A		
Aroclor, 1260	A	A	A		
Aromatic Fuel -50%	A	A	A		
Arsenic Acid	A	A	A	A	A
Arsenic Oxide			A		
Arsenic Tri-chloride	D	C	A		
Arsenic Trioxide	D		A		
Arsenic Trisulfide	D		A		
Arsenites			A		
Arsine			A		
Aryl Orthosilicate			A		
Ascorbic Acid	C	A	A		
Askarel Transformer Oil	A	A	A		
Aspartic Acid	C	B	A		
Asphalt	A	A	A		A
ASTM Oil, No. 1	A	A	A	C	
ASTM Oil, No. 2	A	A	A	D	
ASTM Oil, No. 3	A	A	A	D	
ASTM Oil, No. 4	A	A	A	D	
ASTM Oil, No. 5	A	A	A		
ASTM Reference Fuel A	A	A	A	C	
ASTM Reference Fuel B	A	A	A	B	
ASTM Reference Fuel C	A	A	A	B	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
ASTM Reference Fuel D	A	A	A		
ATL-857	A	A			
Atlantic Dominion F	A	A			
Atlantic Utro Gear-e	A	A			
Atlantic Utro Gear-EP Lube	A	A			
Aure 903R (Mobil)	A	A	A		
AUREX 256		A	A		
Automatic Transmission Fluid	A	A	A	D	
Automotive Brake Fluid	D	B	A		
AXAREL 9100					
Azobenzene			A		
Bardol B	A	A	A		
Barium Carbonate	C	A	A	A	A
Barium Chlorate	C	A	A		
Barium Chloride	A	A	A	A	A
Barium Cyanide	A	A	A		A
Barium Hydroxide	A	A	A	A	A
Barium Iodide	A	A	A		
Barium Nitrate	C	A	A		A
Barium Oxide	A	A	A		
Barium Peroxide	C	B	A		
Barium Polysulfide	C	B	A		
Barium Salts	A	A	A		
Barium Sulfate	A	A	A		A
Barium Sulfide	A	A	A	A	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Bayol 35	A	A	A		
Bayol D	A	A	A		
Beer	A	A	A	A	
Beet Sugar Liquids	A	A	A	A	A
Beet Sugar Liquors	A	A	A	A	
Benzaldehyde	D	C	B	B	A
Benzaldehyde Disulfonic Acid			B		
Benzamide	A	A	A		
Benzanthrone	A	A	A		
Benzene	A	A	A	D	A
Benzene Hexachloride		A	A		
Benzenesulfonic Acid 10%	A	A	B		
Benzidine	A	A	A		
Benzidine 3 Sulfonic Acid	A	A	A		
Benzil	A	A	A		
Benzilic Acid	A	A	A		
Benzine (Ligroin)	A	A	A	D	
Benzocatechol	A	A			
Benzochloride	A	A	A		
Benzoic Acid	A	A	A	B	A
Benzoin	A	A	A	D	
Benzonitrile	C	B	A		
Benzophenone	A	A	A		
Benzoquinone	A	A	A		
Benzotrichloride	A	A	A		
Benzotrifluoride	A	A	A		
Benzoyl Chloride	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Benzoyl Peroxide			A		
Benzoylsulfonilic Acid	A	A	A		
Benzyl Acetate	C	C	A		
Benzyl Alcohol	A	A	A		
Benzyl Amine		C	B		
Benzyl Benzoate	A	A	A		
Benzyl Bromide	A		A		
Benzyl Butyl Phthalate	C	A	A		
Benzyl Chloride	A	A	A		
Benzyl Phenol	A	A	A	D	
Benzyl Salicylate	A		A		
Beryllium Chloride	A		A		
Beryllium Fluoride	A		A		
Beryllium Oxide	A		A		
Beryllium Sulfate	C		A		
Bismuth Carbonate	C	A	A		
Bismuth Nitrate	C		A		
Bismuth Oxide	C		A		
Bittern			A		
Black Liquor	A	A	A		
Black Point 77	A	A			
Blast Furnace Gas	A	A	A		
Bleach Liquor	A	A	A		
Bleach Solutions	A		B		
Borax	A		A		A
Borax Solutions	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Bordeaux Mixture	A	A	A		
Boric Acid	A	A	A	A	A
Boric Oxide	C		A		
Borneol	A	A	A		
Bornyl Acetate	A	B	A		
Bornyl Chloride	A		A		
Bornyl Formate	A		A		
Boron Fluids (HEF)	A	A	A		
Boron Hydride			A		
Boron Phosphate			A		
Boron Tribromide			A		
Boron Trichloride			B		
Boron Trifluoride			B		
Boron Trioxide			A		
Brake Fluid DOT3 (Glycol Type)	D	A	A		
Bray GG-130	A				
Brayco 719-R (VV-H-910)	D	A			
Brayco 885 (MIL-L-6085A)	A	A			
Brayco 910	D	B			
Bret 710	D				
Brine	A	A	A		
Brine (Seawater)	A				
Brom - 113					
Brom - 114	B				
Bromic Acid	C		A		
Bromine	A	A	A		D
Bromine Pentafluoride Factory	D	C	B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Bromine Tri-fluoride Factory	D	C	B		
Bromine Water	A		B		
Bromobenzene	A	A	A		
Bromobenzene Cyanide	C		A		
Bromochlorotri-fluoroethane (Halothane)	A	A	A		
Bromoform	A		A		
Bromomethane (Methyl Bromide)	A	A	A		
Bromotrifluoro-ethylene (BFE)		A			
Bromotrifluoro-methane (F-13B1)		B	B		
Brucine Sulfate	C		A		
Buffered Oxide Etchants					
Bunker Oil	A	A	A		
Bunker's "C" (Fuel Oil)	A	A	A		
Butadiene (Mon-omer)	A	A	A		A
Butane	A	A	A	A	A
Butane, 2, 2-Dimethyl	A		A		
Butane, 2, 3-Dimethyl	A		A		
Butanediol	C	A	B		
Butanol (Butyl Alcohol)	A	A	A		A
Butene 2-Ethyl (1-Butene 2-Ethyl)	A		A		
Butter-Animal Fat	A				A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Butyl Acetate or n-Butyl Acetate	D	A	A		C
Butyl Acetyl Ricinoleate	A		A		
Butyl Acrylate	D	C	A		
Butyl Alcohol	A	A	A		
Butyl Alcohol (Secondary)	A	A	A		
Butyl Alcohol (Tertiary)	A	A	A	D	
Butyl Amine or N-Butyl Amine	D	C	A		
Butyl Benzoate	C	A	A		
Butyl Benzoate or n-Butyl Benzoate	A		A		
Butyl Benzolate			A		
Butyl Butyrate or n-Butyl Butyrate	A	A	A		
Butyl Carbitol	C	A	A		
Butyl Cellosolve	D		A		
Butyl Cellosolve Acetate	C		A		
Butyl Cellosolve Adipate	B		A		
Butyl Chloride	A		A		
Butyl Ether or n-Butyl Ether	D		A		
Butyl Glycolate	C		A		
Butyl Lactate	C		A		
Butyl Laurate	C		A		
Butyl Mercaptan (Tertiary)	A		A		
Butyl Methacrylate	C		A		
Butyl Oleate	A		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Butyl Oxalate	C		A		
Butyl Stearate	A		A		
Butylbenzoic Acid	A		A		
Butylene	A	A	A		A
Butyraldehyde	D		B		
Butyric Acid	B		A		A
Butyric Anhydride	C		A		
Butyrolacetone	C		A		
Butyryl Chloride	A		A		
Cadmium Chloride	C		A		
Cadmium Cy- nide	C	A	A		
Cadmium Nitrate	C		A		
Cadmium Oxide	C		A		
Cadmium Sulfate	C		A		
Cadmium Sulfide	C		A		
Calcine Liquors	A		A		
Calcium Acetate	D	A	A		
Calcium Arsenate	C		A		
Calcium Benzoate	A		A		
Calcium Bicar- bonate	C	A	A		
Calcium Bisulfide	C	A	A		B
Calcium Bisulfite	B		A		A
Calcium Bromide	A		A		
Calcium Carbide		A	A		
Calcium Carbo- nate	A	A	A	A	A
Calcium Chlorate	C	A	A		A
Calcium Chloride	A	A	A	A	D
Calcium Chro- mate	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Calcium Cyanamide			A		
Calcium Cyanide		A	A		
Calcium Fluoride	A		A		
Calcium Gluconate	C		A		
Calcium Hydride	A		A		
Calcium Hydro-sulfide	C		A		
Calcium Hydroxide	A	A	A	A	A
Calcium Hypochlorite	A	A	A	A	C
Calcium Hypophosphite	C		A		
Calcium Lactate	C		A		
Calcium Naphthenate			A		
Calcium Nitrate	A	A	A	A	
Calcium Oxalate	C		A		
Calcium Oxide	A	A	A		
Calcium Permanganate			A		
Calcium Peroxide			A		
Calcium Phenol-sulfonate	C		A		
Calcium Phosphate	A	A	A		
Calcium Phosphate Acid	C		A		
Calcium Propionate	C		A		
Calcium Pyridine Sulfonate			A		
Calcium Salts	A	A	A	A	
Calcium Silicate	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Calcium Stearate	A	A	A		
Calcium Sulfamate	A		A		
Calcium Sulfate	C	A	A	A	A
Calcium Sulfide	A	A	A		
Calcium Sulfite	A	A	A		
Calcium Thiocyanate	C		A		
Calcium Thiosulfate	A	A	A		
Calcium Tungstate	C				
Caliche Liquors	A	A	A		
Camphene	A	A	A		
Camphor	A	A	A		
Camphoric Acid	A	A	A		
Cane Sugar Liquors	A	A	A		A
Capric Acid	A	A	A		
Caproic Acid	A	A	A		
Caproic Aldehyde	D	C	A		
Caprolactam	A	C	A		
Capronaldehyde	A		B		
Carbamate	A	A	A		
Carbazole			A		
Carbitol	B	B	A		
Carbolic Acid (Phenol)	A	A	A		
Carbon Bisulfide	A		A		
Carbon Dioxide	A	A	A		A
Carbon Dioxide (Explosive Decompression Use)	A	A	A		A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Carbon Disulfide	A	A	A	D	A
Carbon Fluorides	A		A		
Carbon Mon-oxide	A	A	A	A	A
Carbon Tetrabromide			A		
Carbon Tetrachloride	A	A	A	D	B
Carbon Tetrafluoride	A		B		
Carbonic Acid	A	A	A	A	B
Casein	C	A	A		
Castor Oil	A	A	A	C	
Caustic Lime	C	A	A		
Caustic Potash	C	A	A		
Caustic Soda (Sodium Hydroxide)	C	B	A		
Cellosolve	D	C	A	C	
Cellosolve, Acetate	D	C	A	C	
Cellosolve, Butyl	D		A		
Celluguard	A	A			
Cellulose Acetate	C	B	A		
Cellulose Acetate Butyrate	C		A		
Cellulose Ether	C		A		
Cellulose Nitrate*	C		A		
Cellulose Tripropionate	C		A		
Cellulube (Phosphate Esters)		A	A		
Cellutherm 2505A	A	A			
Cerium Sulfate	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Cerous Chloride	C		A		
Cerous Fluoride	C		A		
Cerous Nitrate	C		A		
Cetane (Hexadecane)	A	A	A		
Cetyl Alcohol	A	A	A		
Chaulmoogric Acid			A		
China Wood Oil (Tung Oil)	A	A	A		
Chloral	C	C	A		
Chloramine		C	A		
Chloranthraquinone	A		A		
Chlordane	A	A	A		
Chlorextol	A	A	A		
Chloric Acid	C		A		D
Chlorinated Solvents, Dry	A	A	A		
Chlorinated Solvents, Wet	A	A	A		A
Chlorine (Dry)	A	A	A		A
Chlorine (Plasma)		A	B		
Chlorine (Wet)		A	B		A
Chlorine Dioxide	A	A	A		
Chlorine Dioxide, 8% Cl as NaClO ₂ in solution	A	A	B		
Chlorine Tri-fluoride	D	C	B		
Chlorine Water	A	A	A		
Chloro 1-Nitro Ethane (1-Chloro 1-Nitro Ethane)	D	C			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Chloro Oxyfluorides			B		
Chloro Xylenols	A		A		
Chloroacetaldehyde	C	C	B		
Chloroacetic Acid	D	B	A	B	
Chloroacetone	D	B	A		
Chloroacetyl Chloride			A		
Chloroamino Benzoic Acid	C		A		
Chloroaniline	C		A		
Chlorobenzaldehyde	C	C	A		
Chlorobenzene	A	A	B	D	A
Chlorobenzene (Mono)	A		A		
Chlorobenzene Chloride	A		A		
Chlorobenzene Trifluoride	A		A		
Chlorobenzo-chloride	A		A		
Chlorobenzotri-fluoride	A		A		
Chlorobromo Methane	A	A	A		
Chlorobromopropane	A	A	A		
Chlorobutadiene	A	A	A		
Chlorobutane (Butyl Chloride)	A	A	A		
Chlorododecane	A	A	A		
Chloroethane	A	A	A		
Chloroethane Sulfonic Acid	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Chloroethylbenzene	A	A	A		
Chloroform	A	A	A	C	A
Chlorohydrin	C		A		
Chloronaphthalene or o-Chloronaphthalene	A	A	A		
Chloronitrobenzene	C		A		
Chlorophenol or o-Chlorophenol	A	A	A		
Chloropicrin	A	A	A		
Chloroprene	A	A	A		
Chlorosilanes	X		A		
Chlorosulfonic Acid	D		A	D	
Chlorotoluene	A	A	A		
Chlorotoluene Sulfonic Acid	C		A		
Chlorotoluidine	A	B	A		
Chlorotrifluoroethylene (CTFE)			B		
Chlorox	A	A	A		A
Chloroxylols	X				
Cholesterol	A		A		
Chrome Alum	A	A	A		
Chrome Plating Solutions	A	A	A		
Chromic Acid	A	A	A	A	B
Chromic Chloride			A		
Chromic Fluorides			A		
Chromic Hydroxide			A		
Chromic Nitrates			A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Chromic Oxide	A	A	A		
Chromic Phosphate			A		
Chromic Sulfate			A		
Chromium Potassium Sulfate (Alum)	A		A		
Chromyl Chlorides			A		
Cinnamic Acid	A	A	A		
Cinnamic Alcohol	A	A	A		
Cinnamic Aldehyde	A	B	A		
Circo Light Process Oil	A		A		
Citric Acid	A	A	A	A	A
City Service #65 #120 #250	A				
City Service Koolmoter-AP Gear Oil 140-EP Lube	A				
City Service Pace-maker #2	A				
Clorox	A	A	A		
Coal Tar	A	A	A		
Cobalt Chloride	A	A	A		
Cobalt Chloride, 2N	A		A		
Cobaltous Acetate	C		A		
Cobaltous Bromide	A		A		
Cobaltous Linoleate			A		
Cobaltous Naphthenate			A		
Cobaltous Sulfate	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Coconut Oil	A	A	A	C	
Cod Liver Oil	A	A	A		
Codeine	A		A		
Coffee	A		A		
Coke Oven Gas	A	A	A		
Coliche Liquors			A		
Convelex 10					
Coolanol 20 25R 35R 40& 45A (Monsanto)	A		A		
Copper Acetate	D	C	A		
Copper Ammonium Acetate	C		A		
Copper Carbonate	C	A	A		
Copper Chloride	A	A	A	A	D
Copper Cyanide	A	A	A		A
Copper Gluconate	C		A		
Copper Naphthenate			A		
Copper Nitrate	A	A	A	A	A
Copper Oxide	A		A		
Copper Salts	A	A	A		
Copper Sulfate	A	A	A	A	A
Copper Sulfate 10%	A	A	A		B
Copper Sulfate 50%	A	A	A		
Corn Oil	A	A	A	A	
Cottonseed Oil	A	A	A		
Creosote, Coal Tar	A	A	A		
Creosote, Wood	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Cresol (Methyl Phenol)	A	A	A		
Cresols	B	A	A	D	A
Cresylic Acid	A	A	A	B	A
Crotonaldehyde	A	B	A		
Crotonic Acid	A	A	A		
Crude Oil	A	A	A		
Cumaldehyde	A		A		
Cumene	A	A	A		
Cumene Hydro-peroxide			A		
Cupric Sulfate	A	A			
Cutting Oil	A	A	A		
Cyanamide			A		
Cyanides			A		
Cyanogen Chloride			B		
Cyanogen Gas			A		
Cyanohydrin			A		
Cyanuric Chloride			A		
Cyclohexane	A	A	A		A
Cyclohexanol	A	A	A		
Cyclohexanone	D	A	A	D	
Cyclohexene	A	B	A	D	
Cyclohexylamine	A	C	B		
Cyclohexylamine Carbonate			B		
Cyclohexylamine Laurate	A	C	B		
Cyclopentadiene	A	B	A		
Cyclopentane	A	A	A		
Cyclopolyolefins	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Cymene or p-Cymene	A	A	A		
DDT (Dichlorodiphenyltrichloroethane)	A	A	A		
Decalin	A	A	A		
Decane	A		A		
Delco Brake Fluid	A		A		
Denatured Alcohol	A	A	A		
Detergent, Water Solution	A	A	A		A
Developing Fluids (Photo)	A	A	A		
Dexron	A	A	A		
Dextrin	A	A	A		
Dextro Lactic Acid	C		A		
Dextron	A				
Dextrose	C	A	A		
DI Water	B	A	B		
Diacetone	D	B	A		
Diacetone Alcohol	D	B	A		
Dialkyl Sulfates	C		A		
Diallyl Ether		B	A		
Diallyl Phthalate		A	A		
Diamylamine	A	C	B		
Diazinon	B	A	A		
Dibenzyl (sym-Diphenylethane)	A	A	A		
Dibenzyl Ether Factory	D	C	A		
Dibenzyl Sebacate	B	A	A		
Diborane			A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Dibromoethane	A		A		
Dibromoethyl Benzene	A	A	A		
Dibutyl Cellosolve Adipate	C		A		
Dibutyl Ether Factory	C	B	A		
Dibutyl Methylene-dithio Glycolate	A		A		
Dibutyl Phthalate	C	A	A	A	
Dibutyl Sebacate	B	A	A		
Dibutyl Thioglycolate	A		A		
Dibutyl Thiourea	A		A		
Dibutylamine	D	C	A		
Dichloroacetic Acid	A	C	A		
Dichloroaniline	C	A	A		
Dichlorobenzene or o-Dichlorobenzene	A	A	A		
Dichlorobenzene or p-Dichlorobenzene	A	A	A		
Dichlorobutane	A	A	A		
Dichlorobutene	A	A	A		
Dichlorodiphenyl-Dichloroethane (DDD)	A		A		
Dichloroethane	A	A	A		A
Dichloroethylene	A	A	A		
Dichlorohydrin	C		A		
Dichloroisopropyl Ether	C	B	A		
Dichloromethane	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Dichlorophenol	A	A	A		
Dichlorophenoxy-acetic Acid	A		A		
Dichloropropane	A	A	A		
Dichloropropene	A	A	A		
Dichlorosilane			A		
Dicyclohexylamine	D	B	A		
Dicyclohexylammonium Nitrate	C				
Dieldrin	A	B	A		
Diesel Oil	A	A	A	D	A
Di-ester Lubricant MIL-L-7808	A	A	A		
Di-ester Synthetic Lubricants	A	A			
Diethanolamine (DEA)	C	B		A	
Diethyl Benzene	A	A	A		
Diethyl Carbonate	C	B	A		
Diethyl Ether	D	B	A		
Diethyl Phthalate	A		A		
Diethyl Sebacate	B		A		
Diethyl Sulfate	C		A		
Diethylamine	D		B		A
Diethylaniline	C		A		
Diethylene Glycol	A	A	A	A	A
Diethylenetriamine		C	A		
Difluorodibromomethane			A		
Difluoroethane	A	D	A		
Difluoromono-chloroethane	A		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Diglycol Chloroformate	C		A		
Diglycolic Acid	C		A		
Dihydroxydiphenylsulfone	C		A		
Diisobutyl Ketone		C	A		
Diisobutylcarbinol	A		A		
Diisobutylene	A		A		
Diisooctyl Sebacate 3 2	B				
Diisopropyl Ether (DIPE)		A	A		
Diisopropyl Ketone	D	C	A		
Diisopropylbenzene	A		A		
Diisopropylidene Acetone	A	C	A		
Dimethyl Acetamide	C	A	A		
Dimethylaniline (Xylidine)	A	B	A		
Dimethyldisulfide (DMDS)	A		A		
Dimethyl Ether	B	C	A		
Dimethyl Formaldehyde	C	C	A		
Dimethyl Formamide (DMF)	D	A	A	B	
Dimethylhydrazine	C	B	A		
Dimethyl Phenyl Carbinol	A		A		
Dimethyl Phenyl Methanol	A		A		
Dimethyl Phthalate	B	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Dimethyl Sulf-oxide (DMSO)	C	C	A	D	
Dimethyl Terephthalate (DMT)	A	A	A		
Dimethylamine (DMA)	D		B		
Dinitrochloroben- zene	A		A		
Dinitrogen Tetr- oxide			B		
Dinitrotoluene (DNT)	D		A		
Diethyl Phthalate	B	A	A	A	
Diethyl Sebacate	B	A	A		
Diethylamine	A		A		
Dioxane	D	C	A	D	
Dioxolane	D	C	A		
Dipentene	A	A	A		
Diphenyl	A	A	A	A	
Diphenyl Oxides	A	A	A		
Diphenylamine (DPA)	A		A		
Diphenylene Oxide		A	A		
Diphenylpropane	A	A	A		
Disilane		A			
Di-Tert-Butyl Per- oxide			A		
Dodecylbenzene	A	A	A		
Dow Chemical 50-4	D	A			
Dow Chemical ET378					
Dow Chemical ET588	D				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Dow Corning -11	A				
Dow Corning -1208, 4050, 6620, F-60, XF-60	A	A			
Dow Corning -1265 Fluorosilicone Fluid	A				
Dow Corning -200	A	A			
Dow Corning -220	A	A			
Dow Corning -3	A	A			
Dow Corning -33	A	A			
Dow Corning -4	A	A			
Dow Corning -44	A	A			
Dow Corning -5	A	A			
Dow Corning -510	A	A			
Dow Corning -55	A	A			
Dow Corning -550	A	A			
Dow Corning -704	A	A			
Dow Corning -705	A	A			
Dow Corning -710	A	A			
Dow Corning F-61	A	A			
Dow Guard	A	A			
Dowanol P Mix			A		
Dowtherm, 209	D	A	A		
Dowtherm, A	A	A	A		
Dowtherm, E	A	A	A		
Drinking Water	A	A	A	A	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Dry Cleaning Fluids	A	A	A		
DTE 20 Series, Mobil	A	A	A		
DTE named series, Mobil, light-heavy	A	A	A		
Elco 28-EP lubricant	A	A			
Epichlorohydrin	D	C	B		
Epoxy Resins	D	B	A		
Erucic Acid					
Esam-6 Fluid	D	B			
Esso Fuel 208	A	A			
Esso Golden Gasoline	A	A			
Esso Motor Oil	A	A			
Esso Transmission Fluid (Type A)	A	A			
Esso WS2812 (MIL-L-7808A)	A	A			
Esso XP90-EP Lubricant	A	A			
Esstic 42, 43	A	A			
Ethane	A	A	A		A
Ethanol	C	A	A		A
Ethanol Amine	D	C	A		A
Ethers	C	C	A	C	A
Ethoxyethyl Acetate (EGMEEA)	C				A
Ethyl Acetate-Organic Ester	D	B	A	B	
Ethyl Acetoacetate	D	C	A		
Ethyl Acrylate	D	C	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ethyl Alcohol	C	A	A	A	
Ethyl Ammonium Dichloride					
Ethyl Benzene	A	A	A		
Ethyl Benzoate	A	A	A	D	
Ethyl Bromide	A	A	A	D	
Ethyl Cellosolve	D	B	A		
Ethyl Cellulose	D	C	A		
Ethyl Chloride	A	A	A	C	A
Ethyl Chlorocarbonate	A	A	A		
Ethyl Chloroformate	D		A		
Ethyl Ether	D	C	A	C	
Ethyl Formate	A	A	A		
Ethyl Hexanol	A	A	A		
Ethyl Lactate	C				
Ethyl Mercaptan X 2	B	A	A		
Ethyl Nitrite	C				
Ethyl Oxalate	B	A	A		
Ethyl Pentachlorobenzene	A	A	A		
Ethyl Pyridine	A				
Ethyl Silicate	A	A	A		
Ethyl Stearate	A				
Ethyl Sulfate	D		A		D
Ethyl Tertiary Butyl Ether		A			
Ethyl Valerate	A				
Ethylacrylic Acid					
Ethylamine	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ethylcyclopentane	A	A			
Ethylene	B	A	A		
Ethylene Chloride	B	A	A		A
Ethylene Chlorohydrin	A	A	A	A	
Ethylene Cyanohydrin	A				
Ethylene Diamine	D	A	B		
Ethylene Dibromide	A	A	A		
Ethylene Dichloride	A	A	A	C	A
Ethylene Glycol	A	A	A	A	A
Ethylene Hydrochloride	A				
Ethylene Oxide	D	D	A	A	A
Ethylene Oxide, (12%) and Freon 12 (80%)	D				
Ethylene Trichloride	A	A	A		
Ethyleneimine					
Ethylmorpholine Stannous Octoate (50/50 mixture)	D				
Ethylmorpholine	A				
Ethylsulfuric Acid	C				
F-60 Fluid (Dow Corning)	A	A			
F-61 Fluid (Dow Corning)	A	A			
Fatty Acids	A	A	A		
FC-43 Heptacosofluorotributylamine	A				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
FC75 & FC77 (Fluorocarbon)	B				
Ferric Acetate	C	C			
Ferric Ammonium Sulfate	C				
Ferric Chloride	A	A	A	A	
Ferric Ferrocyanide	C				
Ferric Hydroxide	C	A	A		
Ferric Nitrate	A	A	A	A	A
Ferric Persulfate	A				
Ferric Sulfate	A	A	A		A
Ferrous Ammonium Citrate	C				
Ferrous Ammonium Sulfate	C	A			
Ferrous Carbonate	C				
Ferrous Chloride		A	A		D
Ferrous Iodide	C				
Ferrous Sulfate	C	A	A		C
Ferrous Tartrate	C				
Fish Oil	A	A	A		
Fisher Reagent					
Fluorinated Cyclic Ethers					
Fluorine (Gas)		A	B	D	
Fluorine (Liquid)	B	A	B		D
Fluorobenzene	A	A	A		
Fluoroboric Acid					
Fluorocarbon Oils					
Fluoroform (Trifluoromethane)					
Fluorolube	B		B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Fluorophosphoric Acid					
Fluorosilicic Acid	B	A	A		B
Fluorosulfonic Acid					
Formaldehyde	D	C	A	D	A
Formamide	C		A		
Formic Acid	D	C	B		B
Freon, 11	B	B	B	A	A
Freon, 112 (Tetra-chlorodifluoroethane)	A	B	B		
Freon, 113	B	C	C		A
Freon, 113 + High and Low Aniline Oil		C			
Freon, 114	A	B	C		
Freon, 114B2	B	B	C		
Freon, 115, 116	B	B	C		
Freon, 12	C	B	B	A	D
Freon, 12 and ASTM Oil #2 (50/50 Mixture)	A	B			
Freon, 12 and Suniso 4G (50/50 Mixture)	A	B			
Freon, 123 (Dichlorotri-fluoroethane)					
Freon, 124 (Chlorotetra-fluoroethane)					
Freon, 125 (Pentafluoroethane)					
Freon, 13	A		A		
Freon, 134a (Tetrafluoroethane)			B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Freon, 13B1	A		B		
Freon, 14	A		A		
Freon, 141b (Dichlorofluoroethane)					
Freon, 142b	B				
Freon, 152a (Difluoroethane)					
Freon, 21	D	D	B		
Freon, 218	A				
Freon, 22 (Chlorodifluoroethane)	D	D	B	A	A
Freon, 22 and ASTM Oil #2 (50/50 Mixture)	B	D	B		
Freon, 23 (Fluoroform)					
Freon, 31	D	D	B		
Freon, 32	D	D	B		
Freon, 502	B		C		
Freon, BF (R112)	A	B			
Freon, C316	A	B	B		
Freon, C318	B	B	C		
Freon, K-142b	D		C		
Freon, K-152a	D		C		
Freon, MF (R11)	B				
Freon, PCA (R113)	B		C		
Freon, TA	D	D	C		
Freon, TC	B	B	B		
Freon, TF (R113)	B	B			A
Freon, TMC	B	B	B		
Freon, T-P35	B	B	B		
Freon, T-WD602	A	B	B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Fuel Oil, #6	A	A	A	A	A
Fuel Oil, 1, and 2	A				A
Fuel Oil, Acidic	A				A
Fumaric Acid	A	A	A		
Fuming Sulphuric Acid (20/25% Oleum)	A				
Furaldehyde	D	C			
Furan (Furfuran)	C	C	A		A
Furfural (Furfuraldehyde)	D	C	B	D	A
Furfuraldehyde	D	C			
Furfuryl Alcohol	B	A	A		
Furoic Acid					
Furyl Carbinol					
Fyrquel 150 220 300 550	A	A			
Fyrquel 90, 100, 500	A	A			
Fyrquel A60	D	A			
Gallic Acid	A	A	A	B	A
Gasoline	A	A	A		A
Gelatin	A	A	A	A	A
Germane (Germanium Tetrahydride)					
Girling Brake Fluid	D	B			
Glauber's Salt	A	A	A		
Gluconic Acid	C	B	A		
Glucose	A	A	A	A	A
Glue		A			A
Glutamic Acid	C	B			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Glycerine (Glycerol)	A	A	A	A	A
Glycerol Dichlorohydrin	C				
Glycerol Monochlorohydrin	C				
Glycerol Triacetate	C				
Glycerophosphoric Acid	C				
Glyceryl Phosphate	C				
Glycidol	C				
Glycol Monoether					
Glycolic Acid	C		A	B	
Glycols	A	A	A		
Glyoxylic Acid	C				
Grease Petroleum Base	A		A		A
Green Sulfate Liquor	A	A	B		
Gulf Endurance Oils	A	A			
Gulf FR Fluids (Emulsion)	A				
Gulf FR G-Fluids	A				
Gulf FR P-Fluids	B				
Gulf Harmony Oils	A				
Gulf High Temperature Grease	A				
Gulf Legion Oils	A				
Gulf Paramount Oils	A				
Gulf Security Oils	A				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Gulfcrown Grease	A				
Halothane	A	A	B		
Halowax Oil	A	A	B		
Hannifin Lube A	A				
Heavy Water		A	A		
HEF-2 (High Energy Fuel)	A	A	A		
Helium	A	A	A		
Heptachlor	A	A			
Heptachlorobutene	A				
Heptaldehyde (Heptanal)	D	C			
Heptane or n-Heptane	A	A	A	D	A
Heptanoic Acid	A	A			
Hexachloroacetone	C	B			
Hexachlorobutadiene	A		A		
Hexachlorobutene	A				
Hexachloroethane	A	A			
Hexaethyl Tetraphosphate					
Hexafluoroethane (F-116)					
Hexafluoroxylene					
Hexaldehyde or n-Hexaldehyde	D	C			
Hexamethyldisilazane					
Hexamethylene (Cyclohexane)	A				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Hexamethylene Diammonium Adipate	A				
Hexamethylene-diamine	C				
Hexamethylene-tetramine	C	C			
Hexane or n-Hexane	A	A	A	D	A
Hexene-1 or n-Hexene-1	A	A			
Hexone (Methyl Isobutyl Ketone)	C	C	A		
Hexyl Acetate	C	C			
Hexyl Alcohol	A	A	A		
Hexylene Glycol	C	A	A		
Hexylresorcinol	A	A			
High Viscosity Lubricant, H2	A	A			
High Viscosity Lubricant, U4	A	A			
HiLo MS #1	D				
Houghto-Safe 1010 phosphate ester	A	A			
Houghto-Safe 1055 phosphate ester	A	A			
Houghto-Safe 1120 phosphate ester	A	A			
Houghto-Safe 271 (Water & Glycol Base)	B	A			
Houghto-Safe 416 & 500 Series		A			
Houghto-Safe 5040 (Water/Oil emulsion)	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Houghto-Safe 620 Water/Glycol	B	A			
Hydraulic Oil (Petroleum Base, Industrial)	A	A	A		A
Hydraulic Oils (Synthetic Base)	A	A	A		A
Hydrazine	D	C	A	C	A
Hydrazine (Anhydrous)	D	C	B		
Hydrazine Dihydrochloride	C				
Hydrazine Hydrate	C	C			
Hydriodic Acid	A	A			
Hydroabietyl Alcohol					
Hydrobromic Acid	A	A	A		D
Hydrobromic Acid 40%	A				D
Hydrocarbons, Saturated	A				
Hydrochloric Acid (cold) 37%	A	A	A	A	A
Hydrochloric Acid (hot) 37%	A	A	A	A	D
Hydrochloric Acid, 3 Molar to 158°F	A	A	A	A	D
Hydrochloric Acid, Concentrated Room Temp.	A	A	A	A	D
Hydrochloric Acid, Concentrated to 158°F	A	A	A	A	D
Hydrocyanic Acid	A	A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Hydro-Drive MIH-10 (Petro-leum Base)	A				
Hydro-Drive MIH-50 (Petro-leum Base)	A				
Hydrofluoric Acid (Anhydrous)		C	A	D	D
Hydrofluoric Acid (conc.) Cold		C		D	D
Hydrofluoric Acid (conc.) Hot 4 3	C	C		D	D
Hydrofluorosilicic Acid	A	A	A		D
Hydrogen Bro-mide (Anhydrous)			A		
Hydrogen Chloride (Anhy-drous)			A		
Hydrogen Chloride gas	A				
Hydrogen Cya-nide		A			
Hydrogen Fluo-ride			A		
Hydrogen Fluo-ride (Anhydrous)	D	C	A		
Hydrogen Gas, Cold	A		A	A	A
Hydrogen Gas, Hot	A				
Hydrogen Iodide (Anhydrous)					
Hydrogen Per-oxide	A	A	A	A	C
Hydrogen Per-oxide 90%	A	A	A		B
Hydrogen Sele-nide					

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Hydrogen Sulfide, Dry, Cold	D	C	A	A	A
Hydrogen Sulfide, Dry, Hot	D	C			A
Hydrogen Sulfide, Wet, Cold	D	C	A		A
Hydrogen Sulfide, Wet, Hot	D	C			A
Hydrolube-Water/Ethylene Glycol	A	A	A		
Hydrooxycitronellal	A				
Hydroquinol	A				
Hydroquinone	B	A	B		
Hydroxyacetic Acid	C	B	A		
Hydyne	D	D			
Hyjet	D	C			
Hyjet IV and IVA	D	C			
Hyjet S4	D	C			
Hyjet W	D	C			
Hypochlorous Acid	A	A	A		
Indole	A		A		
Industron FF44	A				
Industron FF48	A				
Industron FF53	A				
Industron FF80	A				
Insulin	C				
Iodic Acid	C				
Iodine	A	A			C
Iodine Penta-fluoride	D	C	B		
Iodoform	B	B	A		A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Isoamyl Acetate	C	C	A		
Isoamyl Butyrate	C	B	A		
Isoamyl Valerate	C	A			
Isoboreol	A				
Isobutane	A	A	A		A
Isobutyl Acetate	C	C			
Isobutyl Alcohol	A	A	A	C	A
Isobutyl Chloride	A	A	A		
Isobutyl Ether	D	C			
Isobutyl Methyl Ketone	C	C			
Isobutyl n-Butyrate	A				
Isobutyl Phosphate	C				
Isobutylene	A	A			
Isobutyraldehyde	D				
Isobutyric Acid	D	B			
Isocrotyl Chloride	A				
Isodecanol	A	A			
Isododecane	A	A			
Isoeugenol	A	A			
Isooctane	A	A	A	D	A
Isopentane	A	A	A		
Isophorone (Ketone)	D	A	A		
Isopropanol	A	A	A		
Isopropyl Acetate	D	B	A	B	A
Isopropyl Alcohol	A	A	A	C	A
Isopropyl Chloride	A	A	A		
Isopropyl Ether	D	C	A	C	A
Isopropylacetone	C	B			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Isopropylamine	C	B			
Jet Fuel A	A	A		D	A
JP-10	A	A			
JP-3 (MIL-J-5624)	A	A			A
JP-4 (MIL-T-5624)	A	A			A
JP-5 (MIL-T-5624)	A	A			A
JP-6 (MIL-J-25656)	A	A			A
JP-8 (MIL-T-83133)	A	A			
JP-9 (MIL-F-81912)	A	A			
JP-9 -11	A	A			
JPX (MIL-F-25604)	D	C			
Kel F Liquids	B		B		
Kerosene (Similar to RP-1 and JP-1)	A	A	A		A
Ketones				D	A
Keystone #87HX-Grease	A				
Lacquer Solvents	D	B	A		A
Lacquers	D	B	A		A
Lactams-Amino Acids	D	A	A		
Lactic Acid, Cold	A	A	A	A	A
Lactic Acid, Hot	A	A	A	B	A
Lactones (Cyclic Esters)	D				
Lard Animal Fat	A	A	A	C	A
Lauric Acid	A	A			A
Lavender Oil	A	A			A
LB 135	A				
Lead (Molten)					

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Lead Acetate	D	B	A		A
Lead Arsenate	C				
Lead Azide					
Lead Bromide	C				
Lead Carbonate	C				
Lead Chloride	C	A	A		
Lead Chromate	C		A		
Lead Dioxide	C	A			
Lead Linoleate	C				
Lead Naphthenate					
Lead Nitrate		A	A		A
Lead Oxide	C	A			
Lead Sulfamate	A		A		
Lehigh X1169	A				
Lehigh X1170	A				
Light Grease	A		A		
Ligroin (Petroleum Ether or Benzene)	A	A	A		A
Lime Bleach	A	A	A		
Lime Sulfur	A	A	A		A
Lindol, Hydraulic Fluid (Phosphate ester type)	B	A	A		
Linoleic Acid	B	A	A	C	B
Linseed Oil	A	A	A	C	A
Liquid Oxygen (LOX)	D	D			
Liquid Petroleum Gas (LPG)	A	A	A		
Liquimoly	A		A		
Lithium Bromide (Brine)	C	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Lithium Carbo-nate	C				
Lithium Chloride	C	A	A		A
Lithium Citrate	C				
Lithium Hydroxide	C	B	A		A
Lithium Hypo-chlorite	C				
Lithium Nitrate	C				
Lithium Nitrite	C				
Lithium Perchlo-rate	C				
Lithium Salicylate	C				
Lithopone	C		A		
Lubricating Oils (Crude & Refined)	A	A	A		A
Lubricating Oils (Synthetic base)	A	A	A		A
Lubricating Oils, Di-ester	A	A			A
Lubricating Oils, petroleum base	A	A			A
Lubricating Oils, SAE 10, 20, 30, 40, 50	A	A			A
Lye Solutions	B	A	A		
Magnesium Chloride	A	A	A	A	C
Magnesium Hydroxide	A	A	A	A	A
Magnesium Salts	A	A	A	A	
Magnesium Sul-fite and Sulfate	A	A	A	A	A
Magnesium Trisi-licate	X				
Malathion	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Maleic Acid	A	A	A	C	A
Maleic Anhydride	D	B	A		
Maleic Hydrazide	C				
Malic Acid	A	A	A		A
Mandelic Acid	C	C			
Manganese Acetate	C	B			
Manganese Carbonate	C		A		
Manganese Chloride	C	A	A		
Manganese Dioxide	C				
Manganese Glucuronate	C				
Manganese Hypophosphite	C				
Manganese Linoleate	C				
Manganese Naphthenate	C				
Manganese Phosphate	C				
Manganese Salts		A		A	
Manganese Sulfate	C	A	A		A
Manganous Chloride	C				
Manganous Phosphate	C				
Manganous Sulfate	C	A			
Mannitol	C				
MCS 312	A				
MCS 352	D				
MCS 463	D				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
MDI (Methylene di-p-phenylene isocyanate)	C				
Mercaptan	A	B			
Mercaptobenzo-thiazole (MBT)	A				
Mercuric Acetate	C	B			
Mercuric Chloride	A	A	A		C
Mercuric Cyanide	C	A	A		A
Mercuric Iodide	C				
Mercuric Nitrate	C		A		
Mercuric Sulfate	C	A			
Mercuric Sulfite	C				
Mercurous Nitrate	C	A			
Mercury	A	A	A	A	A
Mercury Chloride	C			A	
Mercury Fulmi-nate	C				
Mercury Salts	C			A	
Mercury Vapors	A				
Mesityl Oxide (Ketone)	D	C	A		
Meta-Cresol	A				
Metaldehyde	C	C			
Meta-Nitroaniline	C				
Meta-Toluidine	A				
Methacrylic Acid	C		A		
Methallyl Chloride	A				
Methane	A	A			A
Methanol	D	A		A	A
Methoxychlor					

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Methoxyethanol (DGMMA)	C				
Methyl Abietate	A	A			
Methyl Acetate	D	B	A	B	A
Methyl Acetoacetate	D	B	A		
Methyl Acetophenone*	A	B			
Methyl Acrylate	D	C	A		
Methyl Alcohol	D	A	A		
Methyl Amylketone	C	B			
Methyl Anthranilate	A	A			
Methyl Benzoate	A	A	A		
Methyl Bromide	A	A	A	C	B
Methyl Butyl Ketone	D	B	A		A
Methyl Butyrate Cellosolve	C		A		
Methyl Butyrate Chloride	C				
Methyl Carbonate	A	A	A		
Methyl Cellosolve	D	B	A		A
Methyl Cellulose	D	B	A		
Methyl Chloride	A	A	A	C	A
Methyl Chloroacetate	C	B			
Methyl Chloroform	A	A			A
Methyl Chloroformate	A	A	A		
Methyl Chlorosilanes					

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Methyl Cyanide (Acetonitrile)	C		A		
Methyl Cyclohexanone	A	C			
Methyl Dichloride	A		A		
Methyl Ether	A	B	A		
Methyl Ethyl Ketone (MEK)	D	B	A	D	A
Methyl Ethyl Ketone Peroxide	D		A		
Methyl Ethyl Oleate	A				
Methyl Formate		B	A		A
Methyl Hexyl Ketone (2-Octanone)	C	B			
Methyl Iodide	A		A		
Methyl Isobutyl Ketone (MIBK)	D	B	A	D	A
Methyl Isocyanate	C	B			
Methyl Isopropyl Ketone	D	B			A
Methyl Isovalerate	A				
Methyl Lactate	C				
Methyl Mercaptan		A	A		
Methyl Methacrylate	D	B	A	D	
Methyl Oleate	A	A	A		
Methyl Pentadiene	A				
Methyl Phenylacetate	A	B			
Methyl Salicylate		A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Methyl Tertiary Butyl Ether (MTBE)	C	B			
Methyl Valerate	A				
N-Methyl-2-Pyrrolidone (NMP)	C		A	A	
Methylacrylic Acid	C				
Methylal					
Methylamine	C				A
Methylamyl Acetate	C				
Methylcyclopentane	A				
Methylene Bromide	A		A		
Methylene Chloride	B		A	C	A
Methylene Iodide	A				
Methylglycerol	C				
Methylisobutyl Carbinol	A				
Methylpyrrolidine	A				
Methylpyrrolidone	A				
Methylsulfuric Acid	C				
MIL-A-6091	A	A			
MIL-C-4339	A	A			
MIL-C-7024	A	A			
MIL-C-8188	B	A			
MIL-E-9500	A	A			
MIL-F-16884	A	A			
MIL-F-17111	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
MIL-F-25558 (RJ-1)	A	A			
MIL-F-25656	A	A			
MIL-F-5566	A	A			
MIL-F-81912 (JP-9)	A	A			
MIL-F-82522 (RJ-4)	A	A			
MIL-G-10924	A	A			
MIL-G-15793	A	A			
MIL-G-21568	A	A			
MIL-G-25013	A	A			
MIL-G-25537	A	A			
MIL-G-25760	A	A			
MIL-G-3278	A	A			
MIL-G-3545	A	A			
MIL-G-4343	A	A			
MIL-G-5572	A	A			
MIL-G-7118	A	A			
MIL-G-7187	A	A			
MIL-G-7421	A	A			
MIL-G-7711	A	A			
MIL-H-13910	A	A			
MIL-H-19457	A	A			
MIL-H-22251		A			
MIL-H-27601	A	A			
MIL-H-46170 -15°F to +400°F	A	A			
MIL-H-46170 -20°F to +275°F	A	A			
MIL-H-46170 -55°F to +275°F	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
MIL-H-46170 -65°F to +275°F	A	A			
MIL-H-5606 -65°F to +235°F	A	A			
MIL-H-5606 -65°F to +275°F	A	A			
MIL-H-6083	A	A			
MIL-H-7083	B	A			
MIL-H-8446 (MLO-8515)	A	A			
MIL-J-5161	A	A			
Milk	A	A			A
MIL-L-15016	A	A			
MIL-L-15017	A	A			
MIL-L-17331	A	A			
MIL-L-2104	A	A			
MIL-L-21260	A	A			
MIL-L-23699	A	A			
MIL-L-25681	A	A			
MIL-L-3150	A	A			
MIL-L-6081	A	A			
MIL-L-6082	A	A			
MIL-L-6085	A	A			
MIL-L-6387	A	A			
MIL-L-7808	A	A			
MIL-L-7870	A	A			
MIL-L-9000	A	A			
MIL-L-9236	A	A			
MIL-O-3503	A	A			
MIL-P-27402		A			
MIL-R-25576 (RP-1)	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
MIL-S-3136, Type I Fuel	A	A			
MIL-S-3136, Type II Fuel	A	A			
MIL-S-3136, Type III Fuel	A	A			
MIL-S-3136, Type IV Oil High Swell	A	A			
MIL-S-3136, Type IV Oil Low Swell	A	A			
MIL-S-3136, Type V Oil Medium Swell	A	A			
MIL-S-81087	A	A			
MIL-T-5624, JP-4, JP-5	A	A			
MIL-T-83133	A	A			
Mineral Oils	A	A	A		A
Mixed Acids	C				
MLO-7277 Hydr.	A	A			
MLO-7557	A	A			
MLO-8200 Hydr.	A	A			
MLO-8515	A	A			
Mobil 24dte	A	A			
Mobil 254 Lubricant		A			
Mobil Delvac 1100, 1110, 1120, 1130	A	A			
Mobil HF	A	A			
Mobil Nivac 20, 30	A	A			
Mobil SHC 500 Series	A	A			
Mobil SHC 600 Series	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Mobil Therm 600	A	A			
Mobil Velocite c	A	A			
Mobilgas WA200 ATF	A	A			
Mobilgear 600 Series	A	A			
Mobilgear SHC ISO Series	A	A			
Mobilgrease HP	A	A			
Mobilgrease HTS	A	A			
Mobilgrease SM	A	A			
Mobilith AW Series	A	A			
Mobilith SHC Series	A	A			
Mobiljet II Lubricant		A			
Mobilmistlube Series	A	A			
Mobiloil SAE 20	A	A			
Mobilux	A				
Molybdenum Disulfide Grease	A	A			
Molybdenum Oxide	C				
Molybdenum Tri-oxide	C				
Molybdic Acid	C				
Monobromobenzene	A	A	A		
Monobromotoluene	A				
Monobutyl Para-cresol					
Monochloroacetic Acid	C	B	A		C

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Monochlorobenzene	A	A	A		A
Monochlorobutene	A				
Monochlorohydrin					
Monoethanolamine (MEA)	D		A	C	A
Monoethyl Amine	C		A		
Monoisopropylamine	C		A		
Monomethyl Aniline	B				
Monomethyl Ether (Dimethyl Ether)			A		
Monomethyl Ether (Methyl Ether)	A	A	A		
Monomethyl Hydrazine			A		
Monomethylamine (MMA)	C				
Monomethylaniline	B	A			
Mononitrotoluene	C				
Mononitrotoluene & Dinitrotoluene (40/60 Mixture)	C				
Monovinyl Acetylene	A	A			
Mopar Brake Fluid	D				
Morpholine	A		A		B
Motor Oils	A			D	A
Mustard Gas		A	A		A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Myristic Acid	A	A			
Naphthalene	A	A	A	D	A
Naphthalene Chloride	A	A			B
Naphthalene Sul-fonic Acid	A				
Naphthalenic Acid	A	A	A		
Naphthalonic Acid	A				
Naphthenic Acid	A				
Naphthylamine	C	A			
Naptha	A			D	A
Natural Gas	A	A	A		A
Neatsfoot Oil	A		A		
Neon	A	A	A		
Neville Acid	A	A	A		
Nickel Acetate	D	B	A		
Nickel Ammo-nium Sulfate	C	A			A
Nickel Chloride	A	A	A	A	A
Nickel Cyanide	C				
Nickel Nitrate	C	A	A	A	A
Nickel Salts	A		A	A	
Nickel Sulfate	A	A	A	A	B
Nicotinamide (Niacinamide)	A				
Nicotinamide Hydrochloride	C				
Nicotine	A	A			A
Nicotine Sulfate	C				
Niter Cake	A	A	A		
Nitric Acid, Red Fuming	B	A	B		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Nitric Acid, White Fuming		A			
Nitric Acid (0 - 50%)	A	A		A	A
Nitric Acid (50 - 100%) 4 3	C	A			B
Nitric Acid 3 Molar to 158°F	C	A	A		
Nitric Acid Concentrated Room Temp.	B	A	A		B
Nitric Acid Concentrated to 158°F	D	A			C
Nitroaniline	C				
Nitrobenzene	B	A	A	D	B
Nitrobenzoic Acid	C				
Nitrocellulose	C				
Nitrochlorobenzene	C				
Nitrochloroform	C				
Nitrodiethylaniline	C				
Nitrodiphenyl Ether					
Nitroethane	D	C	A		
Nitrofluorobenzene	C				
Nitrogen	A	A	A		A
Nitrogen Oxides	C	C			
Nitrogen Tetroxide (N2O4)	D	C	A		
Nitrogen Trifluoride					
Nitroglycerine	C				A
Nitroglycerol	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Nitroisopropylbenzene	C				
Nitromethane	D		A	D	
Nitrophenol	C	C			
Nitropropane	D	C	A		
Nitrosyl Chloride					
Nitrosylsulfuric Acid					
Nitrothiophene	C				
Nitrotoluene	C	B			
Nitrous Acid	C		A	A	B
Nitrous Oxide	A	B			A
Nonane	A	A			
Noryl GE Phenolic					
Nyvac FR200 Mobil	A				
Octachloro Toluene	A	A	A		
Octadecane	A	A	A		
Octanal (n-Octanaldehyde)	A	C	A		
Octane or n-Octane	A	A			A
Octyl Acetate	C	B	A		
Octyl Alcohol	A	A	A		
Octyl Chloride	A	A			
Octyl Phthalate	A	A			
Olefins	A	A			
Oleic Acid	B	A	A	C	A
Oleum (Fuming Sulfuric Acid)	A	A	A	A	A
Oleum Spirits	A	A			
Oleyl Alcohol	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Olive Oil	A	A	A	A	A
Oronite 8200	A	A			
Oronite 8515	A	A			
Ortho-Chloro Ethyl Benzene	A	A	A		
Ortho-Chloroaniline	C				
Ortho-Chlorophenol	C				
Ortho-Cresol	C				
Ortho-Dichlorobenzene	A			D	
Ortho-Nitrotoluene	C				
Orthophos Acid					
OS 45 Type III (OS45)	A	A			
OS 45 Type IV (OS45-1)	A	A			
OS 70	A	A			
Oxalic Acid	A	A	A	B	B
Oxygen, 200°-300°F (Evalute for specific applications)	B	A	A		A
Oxygen, 300°-400°F (Evalute for specific applications)	B	A	A		A
Oxygen, Cold (Evalute for specific applications)	A	A		A	A
Oxygen, Liquid	D	D	D		A
Ozonated Deionized Water	C				
Ozone	A	A	A	A	A
Paint Thinner, Duco	B	A	A	B	

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Palmitic Acid	A	A	A	C	A
Para-Aminobenzoic Acid	C				
Para-Bromobenzylphenyl Ether					
Para-Chlorophenol	C				
Paracymene	A				
Para-Dichlorobenzene	A		A		
Paraffins	A	A	A	D	A
Para-Formaldehyde	C				
Paraldehyde	C	C	A		
Par-al-Ketone	D	C	A		
Para-Nitroaniline	C				
Para-Nitrobenzoic Acid	C				
Para-Nitrophenol	C				
Parathion	A				
Para-Toluene Sulfonic Acid	C				
Parker O Lube	A				
Peanut Oil	A	A	A		A
Pectin (Liquor)	A				
Pelagonic Acid	A				
Penicillin (Liquid)	A		A		
Pentachloroethane	A	A	A		
Pentachlorophenol	C	A	A		
Pentaerythritol	C				
Pentaerythritol Tetranitrate	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Pentafluoroethane (F-125)					
Pentane or n-Pentane	A	A	A		A
Pentane, 2 Methyl	A				
Pentane, 2-4 dimethyl	A				
Pentane, 3-Methyl	A				
Pentoxone		C			
Pentyl Penta-noate	A				
Peracetic Acid	C				
Perchloric Acid - 2N	A	A	A	A	A
Perchloroethylene	A	A	A	C	A
Perfluoropropane					
Perfluorotriethyl-amine					
Permanganic Acid					
Persulfuric Acid (Caro's Acid)					
Petrolatum	A	A	A		A
Petrolatum Ether	A	A			
Petroleum Oil, Above 250°F	B	A	A		
Petroleum Oil, Below 250°F	A	A	A		
Petroleum Oil, Crude	A	A	A		
Phenol	A	A	A	A	A
Phenol, 70% / 30% H2O	A	A	A	A	A
Phenol, 85% / 15% H2O	A	A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Phenolic Sulfo-nate	C				
Phenolsulfonic Acid	C	A	A		
Phenylacetamide	A				
Phenylacetate	C	B	A		
Phenylacetic Acid	C				
Phenylbenzene	A	A	A		
Phenylene Dia-mine		C			
Phenylethyl Alcohol	A				
Phenylethyl Ether	D	B	A		
Phenylethyl Malonic Ester*	A				
Phenylglycerine	C				
Phenylhydrazine	A	B	A		
Phenylhydrazine Hydrochloride	C				
Phenylmercuric Acetate	C	C			
Phorone	D	C	A		
Phosgene		D	A		
Phosphine					
Phosphoric Acid 3 Molar to 158°F	A	A	A	B	
Phosphoric Acid Concentrated Room Temp	A	A	A	B	C
Phosphoric Acid Concentrated to 158°F	A	A	A	B	C
Phosphoric Acid, 20%		A	A	B	A
Phosphoric Acid, 45%		A	A	B	B

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Phosphorus (Molten)					
Phosphorus Oxy-chloride			A		
Phosphorus Tri-chloride	A	A	A		
Phosphorus Tri-chloride Acid	A		A		A
Phthalic Acid	C	A		A	A
Phthalic Anhydride	C	C	A	A	B
Pickling Solution 3 2	B	A			
Picric Acid (aq)	A	A	A	D	A
Picric Acid Molten	A				
Pine Oil	A	A	A		A
Pine Tar	A	A			
Pinene	A	A	A		
Piperazine	A				
Piperidine	A		A		
Piranha (H ₂ SO ₄ :H ₂ O ₂) (70:30)					
Plating Solution (Co,Cu,Au,In,Fe,Pb,Ni,Ag,Sn,Zn)	A	A			A
Plating Solutions Chrome	A	A	A		A
Plating Solutions Others	A	A			A
Pneumatic Service	A		A		
Polyethylene Glycol	C	A			
Polyglycerol	C				
Polyglycol	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Polyvinyl Acetate Emulsion					B
Potassium (Molten)					
Potassium Acetate	D	C	A		
Potassium Acid Sulfate	C				
Potassium Alum	C				
Potassium Aluminum Sulfate	C				A
Potassium Antimonate	C				
Potassium Bicarbonate	C	A	A		A
Potassium Bichromate	C				A
Potassium Bifluoride	C				
Potassium Bisulfate	C	A	A		A
Potassium Bisulfite	C	A	A		
Potassium Bitartrate	C				
Potassium Bromide	C	A	A		A
Potassium Carbonate	C	A	A		A
Potassium Chlorate	C	A	A		A
Potassium Chloride	A	A	A		A
Potassium Chromates	C	A	A		B
Potassium Citrate	C				
Potassium Cupro Cyanide	A				A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Potassium Cyanate	C		A		B
Potassium Cyanide	A	A			
Potassium Dichromate	A	A	A		A
Potassium Diphosphate	C				
Potassium Ferri-cyanide	C	A	A		
Potassium Fluoride	C	A	A		
Potassium Glucocyanate	C				
Potassium Hydroxide 50%	D	A	A		B
Potassium Hypochlorite	C	B	A		
Potassium Iodate	C		A		
Potassium Iodide	C	A	A		
Potassium Metabisulfate	C				
Potassium Metachromate	C				
Potassium Metasilicate					
Potassium Monochromate	C				
Potassium Nitrate	A	A	A		B
Potassium Nitrite	C		A		
Potassium Oxalate	C		A		
Potassium Perchlorate	C	A	A		
Potassium Perfluoro Acetate					
Potassium Permanganate	C	A	A		B

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Potassium Peroxide					
Potassium Persulfate	C		A		
Potassium Phosphate (Acid)	C		A		
Potassium Phosphate (Alkaline)	C				
Potassium Phosphate (Di/Tri Basic)	C				
Potassium Pyrosulfate	C				
Potassium Salts	A		A		
Potassium Silicate		A	A		
Potassium Sodium Tartrate	C				
Potassium Stannate	C				
Potassium Stearate	C	A			
Potassium Sulfate	A	A	A		B
Potassium Sulfide	C	A	A		B
Potassium Sulfite	A	A	A		A
Potassium Tartrate	C		A		
Potassium Thiocyanate	C		A		
Potassium Thiosulfate	C				
Potassium Triphosphate	C				A
Prestone Anti-freeze	A	A			
PRL-High Temp. Hydr. Oil	A	A			
Producer Gas	A	A	A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Propane	A	A	A	A	A
Propanol			A		
Propionaldehyde	C	C	A		
Propionic Acid	C	B	A		A
Propionitrile	A	C	A		
Propyl Acetate	D	C	A		A
Propyl Acetone or n-Propyl Acetone	D	B			
Propyl Alcohol	A	A	A		A
Propyl Nitrate	D		A		
Propyl Propionate	C				
Propylamine	C	C			
Propylbenzene	A	A			
Propylene	A	A	A		
Propylene Chloride	A	A	A		
Propylene Chlorohydrin	A				
Propylene Dichloride	A	A	A		
Propylene Glycol	C	A	A	A	A
Propylene Imine	A				
Propylene Oxide	D	D	A	A	A
Pydraul 90e	A	A			
Pydraul, 10E	D	A			
Pydraul, 115E	A	A			
Pydraul, 230C, 312C, 540C, A200	A	A			
Pydraul, 29ELT 30E, 50E, 65E	A	A			
Pyranol Transformer Oil	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Pyridine	A	A	A	C	B
Pyridine Oil	D	C			
Pyridine Sulfate	C				
Pyridine Sulfonic Acid	C				
Pyrogallol (Pyro-gallic Acid)	A		A		A
Pyrogard 42, 43, 55	A	A			
Pyrogard 53, Mobil Phosphate Ester	A	A			
Pyrogard D, Mobil Water-in-Oil Emulsion	D	A			
Pyroligneous Acid	D	C			
Pyrolube	A	A			
Pyrosulfuric Acid	C				
Pyrosulfuryl Chloride	A				
Pyrrole	D	B	A		
Pyruvic Acid	C				
Quinidine	A				
Quinine	A	A	A		
Quinine Bisulfate	C				
Quinine Hydrochloride	C				
Quinine Sulfate	C				
Quinine Tartrate	C				
Quinizarin	A	A			
Quinoline	A	A			
Quinone	A	A	A		A
Radiation (Gamma, 1.0 E +07 Rads)	D	C			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Raffinate	A	A			
Rapeseed Oil	A	A	A		
Red Line 100 Oil	A	A			
Red Oil (MIL-H-5606)	A	A	A		
Resorcinol	C	A			
Rhodium					
Riboflavin	A	A			
Ricinoleic Acid	A	A			
RJ-1 (MIL-F-25558)	A	A	A		
RJ-4 (MIL-F-82522)	A	A			
Rosin	A	A	A		A
RP-1 (MIL-R-25576)	A	A			
Saccharin Solution	C				
Sal Ammoniac	A	A	A		
Salicylic Acid	A	A	A	A	A
Santo Safe 300	A	A			
Sea (Salt) Water	A	A	A		C
Sebacic Acid	C	A			
Selenic Acid	C				
Selenous Acid	C				
Sewage	A	A	A		
SF 1154 GE Silicone Fluid	A	A			
SF1147 GE Silicone Fluid	A	A			
SF96 GE Silicone Fluid	A	A			
Shell 3XF Mine Fluid (Fire resist hydr.)	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Shell Alvania Grease #2	A	A			
Shell Carnea 19 and 29	A	A			
Shell Diala	A	A			
Shell Iru 905	A	AA			
Shell Lo Hydrax 27 and 29	A	A			
Shell Macome 72	A	A			
Shell Tellus #32 (Petroleum Base)	A	A			
Shell Tellus #68	A	A			
Shell Tellus 27 (Petroleum Base)	A	A			
Shell Tellus 33	A	A			
Shell UMF (5% Aromatic)	A	A			
Shellac	C	A			
Silane					
Silicate Esters	A	A	A		
Silicon Fluoride					
Silicon Tetra-chloride			B		
Silicon Tetra-fluoride					
Silicone Greases	A	A	A		
Silicone Oils	A	A	A	C	A
Silver Bromide	C		A		
Silver Chloride	C		A		C
Silver Cyanide	C	A	A		A
Silver Nitrate	A	A	A	A	A
Silver Sulfate	C	A			
Sinclair Opaline CX-EP Lube	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Skelly, Solvent B, C, E	A				
Skydrol 500 B4	D	C	A	D	
Skydrol 7000	B	A	A		
Skydrol LD-4	D	C			
Soap Solutions	A	A		A	A
Socony Mobile Type A	A	A			
Socony Vacuum AMV AC781 (Grease)	A	A			
Socony Vacuum PD959B	A	A			
Soda Ash	A	A			A
Sodium (Molten)		D			
Sodium Acetate	D	A	A	A	A
Sodium Acid Bisulfate	C				
Sodium Acid Fluoride	C				
Sodium Acid Sulfate	C				
Sodium Aluminate	C	A	A		A
Sodium Aluminate Sulfate	C				
Sodium Anthraquinone Disulfate	C				
Sodium Antimonate	C				
Sodium Arsenate	C	A			
Sodium Arsenite	C		A		
Sodium Benzoate	C	A		A	
Sodium Bicarbonate (Baking Soda)	A	A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sodium Bichromate	C		A		A
Sodium Bifluoride	C				
Sodium Bisulfate or Bisulfite	A	A	A		A
Sodium Bisulfide	C				
Sodium Bitartrate	C				
Sodium Borate	A	A	A		A
Sodium Bromate	C		A		
Sodium Bromide	C	A	A		A
Sodium Carbonate (Soda Ash)	A	A	A	A	B
Sodium Chlorate	C	A	A	A	A
Sodium Chloride	A	A	A	A	C
Sodium Chlorite	C		A		
Sodium Chloroacetate	C	B			
Sodium Chromate	C	A	A		A
Sodium Citrate	C		A		
Sodium Cyanamide	C				
Sodium Cyanate	C				
Sodium Cyanide		A	A	A	A
Sodium Diacetate	C	B			
Sodium Diphenyl Sulfonate	C				
Sodium Diphosphate	C				
Sodium Disilicate	C				
Sodium Ethylate	C		A		
Sodium Ferricyanide	C	A	A		A
Sodium Fluoride	C	A	A	A	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sodium Fluorosilicate	C				
Sodium Glutamate	C				
Sodium Hydride					
Sodium Hydrogen Sulfate	C		A		
Sodium Hydrosulfide	C	A			
Sodium Hydrosulfite	C		A		
Sodium Hydroxide, 3 Molar	B	A	A	A	B
Sodium Hypochlorite	A	A	A	A	A
Sodium Hypophosphate	C				
Sodium Hypophosphite	C				
Sodium Hyposulfite	C		A		A
Sodium Iodide	C	A	A		
Sodium Lactate	C		A		
Sodium Metaphosphate	A	A	A		A
Sodium Metasilicate	C	A	A		A
Sodium Methylate	C				
Sodium Monophosphate	C				
Sodium Nitrate		A	A	A	A
Sodium Oleate	C		A		
Sodium Orthosilicate	C				
Sodium Oxalate	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sodium Perborate	A	A	A		C
Sodium Percarbonate	C				
Sodium Perchlorate	C		A		
Sodium Peroxide	A	A	A		A
Sodium Persulfate	C		A		
Sodium Phenolate	C				
Sodium Phenoxide	C				
Sodium Phosphate (Dibasic)	A	A	A		A
Sodium Phosphate (Mono)	A	A	A		A
Sodium Phosphate (Tribasic)	A	A	A		A
Sodium Plumbite	C				
Sodium Pyrophosphate	C		A		
Sodium Resinate	C				
Sodium Salicylate	C				
Sodium Salts	A	A	A	A	
Sodium Sesquioxide	X	A			
Sodium Silicate	A	A	A		B
Sodium Silicofluoride					
Sodium Stannate	C				
Sodium Sulfate	A	A	A	A	A
Sodium Sulfide and Sulfite	A	A	A	A	B
Sodium Sulfocyanide	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sodium Tartrate	C		A		
Sodium Tetraborate	C		A		A
Sodium Tetraphosphate	C				
Sodium Tetrasulfide	C				
Sodium Thioarsenate	C				
Sodium Thiocyanate	C	A			
Sodium Thiosulfate	A	A	A		A
Sodium Trichloroacetate	C				
Sodium Triphosphate	C				
Solvesso 100, 150		A			
Sorbitol	C	A			
Sour Crude Oil	A	B	B		A
Sour Natural Gas	A	B	A		
Sovasol No. 1, 2, and 3	A				
Sovasol No. 73 and 74	A				
Soybean Oil	A	A	A		A
Spry	A				
SR-10 Fuel	A	A			
SR-6 Fuel	A	A			
Standard Oil Mobilube GX90-EP Lube	A				
Stannic Ammonium Chloride	C				
Stannic Chloride	A	A	A		D

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Stannic Chloride, 50%	A		A		
Stannic Tetra-chloride	C				
Stannous Bisul-fate	C				
Stannous Bro-mide	C				
Stannous Chloride (15%)	A	A			C
Stannous Fluo-ride	C	A			
Stannous Sulfate	C				
Stauffer 7700	A	A	A		
Steam Below 400°F	D	A	A		A
Steam, 400° - 500°F	D	C	A		A
Steam, Above 500°F		C	B		A
Stearic Acid		A	A	C	A
Stoddard Solvent	A	A	A		A
Strontium Ace-tate	C	C			
Strontium Carbo-nate	C				
Strontium Chloride	C				
Strontium Hydroxide	C				
Strontium Nitrate	C				
Styrene (Mon-omer)	B	A	A	D	A
Succinic Acid	C	A	A		A
Sucrose Solutions	A	A	A		
Sulfamic Acid	C	A	A		A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sulfanilic Acid	C				
Sulfanilic Chloride	A				
Sulfanilimide	A				
Sulfite Liquors	C	A			A
Sulfolane	B	A			
Sulfonated Oils	A				
Sulfonic Acid	C				
Sulfonyl Chloride	C				
Sulfur	A		A		A
Sulfur (Molten)	A	A			
Sulfur Chloride	A	A	A	D	D
Sulfur Dioxide, Dry	D	A	A		A
Sulfur Dioxide, Liquidified under pressure	D	A	A		A
Sulfur Dioxide, Wet	D	A	A		A
Sulfur Hexa-fluoride	C	A	B		
Sulfur Liquors	C				C
Sulfur Mono-chloride	A				
Sulfur Tetra-fluoride		A			
Sulfur Trioxide Dry	A	A	A		C
Sulfuric Acid (20% Oleum)	C	A	A		C
Sulfuric Acid, 3 Molar to 158°F	A	A	A		D
Sulfuric Acid, Concentrated Room Temp	A	A	A		D

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Sulfuric Acid, Concentrated to 158°F	A	A	A		D
Sulfuric Chlorohydrin (Chlorosulfonic Acid)	C	B			
Sulfurous Acid	A	B	A	A	B
Sunoco #3661	A	A			
Sunoco All purpose grease	A	A			
Sunoco SAE 10	A	A			
Sunsafe (Fire resist. hydr. fluid)	A	A			
Super Shell Gas	A	A			
Surfuryl Chloride	C				
Swan Finch EP Lube	A	A			
Swan Finch Hypoid-90	A				
Tallow	A	A			A
Tannic Acid (10%)	A	A	A		A
Tar, bituminous	A	A	A		A
Tartaric Acid	A	A	A	A	B
Tellone II					
Terephthalic Acid	C	A			
Terpineol	A	A	A		
Terpinyl Acetate	B	B			
Tertiary Amyl Methyl Ether (TAME)					
Tertiary Butyl Catechol or p-tert-butylcatechol	A		A		
Tertiary Butyl Mercaptan	A		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Tetrabromoethane	A	A	A		
Tetrabromomethane	A	A	A		
Tetrabutyl Titanate	A	A	A		
Tetrachloroethylene	A	A	A		
Tetrachloroethane	A	A	A		A
Tetraethyl Lead	A	A	A		
Tetraethyl Lead "Blend" V1164-75	A				
Tetraethyl Orthosilicate (TEOS)					
Tetrahydrofuran	D	B	A	D	A
Tetralin	A	A	A		
Tetramethyl Ammonium Hydroxide	C				
Tetramethylcyclotetrasiloxane (TMCTS)					
Tetramethyldihydropyridine	A				
Tetramethyldihydropyridine	A				
Tetraphosphogluconate	C				
Tetraphosphoric Acid					
Texaco 3450 Gear Oil	A	A			
Texaco Capella A and AA	A	A			
Texaco Meropa 220 (No Lead)	A	A			
Texaco Regal B	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Texaco Uni-Temp Grease	A	A			
Texamatic "A" 1581 Fluid	A	A			
Texamatic "A" 3401 Fluid	A	A			
Texamatic "A" 3525 Fluid	A	A			
Texamatic "A" 3528 Fluid	A	A			
Texamatic "A" Transmission Oil	A	A			
Texas 1500 Oil	A	A			
Therminol 44	A	A			
Therminol 55	A				
Therminol 66					
Therminol FR					
Therminol VP-1, 60, 65	A				
Thio Acid Chloride					
Thioamyl Alcohol	A				
Thiodiacetic Acid	C				
Thioethanol	C				
Thioglycolic Acid	C		A		
Thiokol TP-90B	A	A			
Thiokol TP-95	A	A			
Thionyl Chloride	A	A	A	A	
Thiophene (Thio-furan)	C	C	A		
Thiophosphoryl Chloride	C				
Thiourea	C				
Thorium Nitrate	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Tidewater Multi-gear, 140 EP Lube	A				
Tidewater Oil-Beedol	A				
Tin Ammonium Chloride	C				
Tin Chloride	A			A	
Tin Tetrachloride	A			A	
Titanic Acid	C				
Titanium Dioxide	C	A			
Titanium Sulfate	C		A		
Titanium Tetra-chloride	A	A	B		C
Toluene	A	A	A		A
Toluene Bisodium Sulfite					
Toluene Diisocyanate (TDI)	D	C			
Toluene Sulfonyl Chloride	A				
Toluenesulfonic Acid	C				
Toluidine	A	B	A		
Toluol	C	C	A	D	A
Toluquinone	A				
Tolylaldehyde	C				
Transformer Oil	A	A	A		A
Transmission Fluid Type A	A	A	A		
Triacetin	D	B	A		
Triaryl Phosphate	A				
Tribromomethylbenzene	A				
Tributoxyethyl Phosphate	A		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Tributyl Citrate	C	C			
Tributyl Mercaptan	A	A	A		
Tributyl Phosphate	D	C	A		A
Tributylamine					
Trichloroacetic Acid	C	B	A	B	B
Trichloroacetyl Chloride	A				
Trichlorobenzene	A	A	A		
Trichloroethane	A	A	A	C	A
Trichloroethanolamine	C				
Trichloroethylene	A	A	A	C	A
Trichloromethane	A				
Trichloronitromethane (Chloropicrin)	C				
Trichlorophenylsilane					
Trichloropropane	A	A	A	C	A
Trichlorosilane	A				
Tricresyl Phosphate	B	A		A	A
Triethanol Amine	D	C		C	A
Triethyl Phosphate	A				
Triethylaluminum		A			
Triethylborane		A	A		
Triethylene Glycol	C	A			
Triethylenetetramine	C	C			
Trifluoroacetic Acid	C		A		

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Trifluoroethane	A		A		
Trifluoromethane	A				
Trifluorovinylchloride	A				
Triisopropylbenzylchloride	A				
Trimethylamine	C	C			
Trimethylamine (TMA)	C				
Trimethylbenzene	A				
Trimethylborate (TMB)	A				
Trimethylpentane	A	A	A		
Trinitrololuene (TNT)	B	A	A		
Trioctyl Phosphate	B	A	A		
Triphenylphosphite	C		A		
Tripoly Phosphate	B		A		
Tripotassium Phosphate	C				
Trisodium Phosphate	C				A
Tritium					
Tung Oil (China Wood Oil)	A	A			
Tungsten Hexafluoride					
Tungstic Acid					
Turbine Oil	A	A	A		
Turbine Oil #15 (MIL-L-7808A)	A	A			
Turbo Oil #35	A	A			
Turpentine	A	A	A	D	A

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Type I Fuel (MIL-S-3136)(ASTM Ref. Fuel A)	A	A	A		
Type II Fuel MIL-S-3136	A	A	A		
Type III Fuel MIL-S-3136(ASTM Ref. Fuel B)	A	A	A		
Ucon Hydrolube J-4	A	A			
Ucon Lubricant 50-HB-100	A	A			
Ucon Lubricant 50-HB-260	A	A			
Ucon Lubricant 50-HB-5100	A	A			
Ucon Lubricant 50-HB55	A	A			
Ucon Lubricant 50-HB-660	A	A			
Ucon Lubricant LB-1145	A	A			
Ucon Lubricant LB-135	A	A			
Ucon Lubricant LB-285	A	A			
Ucon Lubricant LB-300X	A	A			
Ucon Lubricant LB-625	A	A			
Ucon Lubricant LB-65	A	A			
Ucon Oil 50-HB-280X	C	C			
Ucon Oil Heat Transfer Fluid 500 (Polyalkalene Glycol)	A	A			
Ucon Oil LB-385	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Ucon Oil LB-400X	A	A			
Undecylenic Acid	A				
Undecylic Acid	A				
Univis 40 (Hydr. Fluid)	A	A			
Univolt #35 (Mineral Oil)	A	A			
Unsymmetrical Dimethyl Hydrazine (UDMH)	D		B		
UPDI(Ultrapure Deionized Water)	C				
Uranium Hexachloride	A				
Uranium Hexafluoride					
Uranium Sulfate Factory					
Uric Acid	C		A	A	
Valeraldehyde	C	C			
Valeric Acid	C	A	A		
Vanadium Oxide	A				
Vanadium Pentoxide	A				
Varnish	A	A	A		A
Vegetable Oil	A	A	A		A
Versilube F44, F55	A	A			
Versilube F-50	A	A			
Vinegar	C	A	A	A	A
Vinyl Acetate	C	C	A	B	
Vinyl Benzene	A		A		
Vinyl Benzoate	A				
Vinyl Chloride	A	A	A		
Vinyl Fluoride	A	A			

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Vinylidene Chloride	A	A			
Vinylpyridine	A				
Vitriol (White)	C				
VV-H-910	A	A			
Wagner 21B Brake Fluid	D	B			
Water	A	A	A	A	A
Wemco C	A	A			
Whiskey and Wines	A	A	A		A
White Liquor	A	A			A
White Oil	A	A	A		A
White Pine Oil	A	A	A		
Wolmar Salt	A	A			
Wood Alcohol	D	A	A		
Wood Oil	A		A		
Xenon	A	A	A		
Xylene	A	A	A		A
Xylidenes-Mixed-Aromatic Amines	D	C	A		
Xylol	A	A		D	
Zeolites	A	A	A		
Zinc Acetate	D	B	A		A
Zinc Ammonium Chloride	C	B	A		
Zinc Chloride	A	A	A	A	B
Zinc Chromate	C	B			
Zinc Cyanide	C	B			
Zinc Diethyldithiocarbamate	C	B			
Zinc Dihydrogen Phosphate	C				

List of chemical resistance

Fluid	FKM / FPM / Viton	Viton® Extreme ETP-S	FFKM / Kalrez®	TPE (LEZ-SAN Santoprene®)	Stainless Steel AISI 316L
Zinc Fluorosilicate					
Zinc Hydrosulfite	C	B	A		A
Zinc Naphthenate					
Zinc Nitrate	A	A	A		A
Zinc Oxide	A	A			
Zinc Phenolsulfonate	C	B			
Zinc Phosphate	A	A	A		
Zinc Salts	A	A	A	A	
Zinc Silicofluoride					
Zinc Stearate	C	A			
Zinc Sulfate	A	A	A		A
Zinc Sulfide	C	B			
Zirconium Nitrate	A	A			

3.4 Table 2

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Acetaldehyde	B	D	A	A
Acetate LMW	—	—	A	A
Acetic acid <5%	A	—	A	B
Acetic acid >5%	A	B	A	B
Acetic anhydride	C	D	A	B
Acetone	C	D	A	A
Acetonitrile	—	D	A	A
Acetyl bromide	—	—	A	—
Acetyl chloride	C	A	A	A
Air	A	A	A*	A

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Aliphatic hydrocarbons	—	—	—	B
Aluminum chloride	B	A	A	D
Aluminum sulfate	A	A	A	B
Alums	A	A	A	—
Ammonia, gas / liquid	C	D	A	B
Ammonium acetate	—	D	A	B
Ammonium carbonate	C	A	A	B
Ammonium chloride	C	A	A	C
Ammonium hydroxide	A	B	A	A
Ammonium nitrate	C	A	A	A
Ammonium phosphate	A	A	A	B
Ammonium sulfate	A	A	A	B
Amyl acetate	D	D	A	A
Amyl alcohol	D	A	A	A
Amyl chloride	D	A	A	A
Aniline	D	B	A	A
Aniline hydrochloride	D	B	A	D
Aqua regia (80 % HCl, 20 % H)	D	B	A	D
Aromatic hydrocarbons	—	A	—	B
Arsenic salts	—	D	—	—
Barium salts	A	A	A	B
Benzaldehyde	B	D	A	B
Benzenesulfonic acid	D	A	A	B
Bleaching liquors	B	A	A	—
Boric acid	A	A	A	B

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Bromine	D	A	A	D
Butane	D	A	A	A
Butanol (butyl alcohol)	B	A	A	A
Butyl acetate	D	D	A	B
Butyric acid	D	B	A	B
Calcium oxide	A	A	A	A
Calcium salts	B	A	A	B
Carbon bisulfide	D	—	A	A
Carbon dioxide	B	A	A*	A
Carbon tetrachloride	D	A	A	B
Chlorine, dry	D	A	A*	A
Chlorine, wet	D	B	A	C
Chloroacetic acid	—	D	A	B
Chlorobenzene	D	A	B	A
Chlorobromomethane	D	A	A	—
Chloroform	D	A	A	A
Chlorosulfonic acid	D	D	A	D
Chromic acid, 30 %	C	A	A	B
Chromium salts	—	—	—	—
Copper salts	A	A	A	B
Cresol	D	A	A	A
Cyclohexane	D	A	A	A
Cyclohexanone	D	D	A	A
Diacetone alcohol	B	D	A	B
Dimethyl formamide	B	D	A	A
Essential oils	C	—	—	—
Ethanol (ethyl alcohol)	A	A	A	A

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Ether	D	D	A	A
Ethyl acetate	B	D	A	B
Ethyl bromide	D	A	A	—
Ethyl chloride	D	A	A	A
Ethylamine	C	D	—	—
Ethylene chlorohydrin	C	A	A	B
Ethylene dichloride	D	A	A	B
Ethylene glycol	A	A	A	B
Ethylene oxide	D	D	A	B
Fatty acids	C	A	A	B
Ferric chloride	B	A	A	D
Ferric sulfate	B	A	A	B
Ferrous chloride	C	A	A	D
Ferrous sulfate	C	A	A	B
Fluoboric acid	A	—	A	B
Fluoroborate salts	—	—	—	—
Fluosilicic acid	D	A	A	C
Formaldehyde	B	D	A	C
Formic acid, 25 %	B	D	A	B
Freon® TMS	—	—	A	—
Gasoline, high-aromatic	D	A	B	A
Gasoline, nonaromatic	D	A	A	—
Glucose	A	A	A	A
Glue, P.V.A.	A	A	A	A
Glycerin	A	A	A	A
Hydriodic acid	—	A	—	—

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Hydrobromic acid, 30 %	D	A	A	D
Hydrochloric acid (dil)	D	A	A	D
Hydrochloric acid (med)	D	A	A	D
Hydrochloric acid (conc)	D	A	A	D
Hydrocyanic acid	C	A	A	B
Hydrocyanic acid, gas, 10 %	C	A	A	—
Hydrofluoric acid, 50 %	D	D	A	D
Hydrofluoric acid, 75 %	D	D	A	D
Hydrogen peroxide (dil)	A	A	A	B
Hydrogen peroxide, 90 %	B	A	A	B
Hypochlorous acid	D	A	A	—
Iodine solutions	C	A	A	D
Iodoform	—	C	—	A
Kerosene	D	A	A	A
Ketones	—	—	A	A
Lacquer solvents	D	D	A	A
Lactic acid, 3-10%	A	A	A	B
Lead acetate	D	D	A	B
Linseed oil	A	A	A	A
Lithium hydroxide	D	C	A	B
Magnesium chloride	A	A	A	B
Magnesium sulfate	A	A	A	B
Malic acid	B	A	A	A
Manganese salts	B	A	A	D

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Mercury salts	—	A	A	B
Methane	D	A	A*	A
Methanol (methyl alcohol)	A	B	A	A
Methyl chloride	D	B	A	A
Methyl ethyl ketone (MEK)	D	D	A	A
Mixed acid (40% H ₂ SO ₄ , 15 % HNO ₃)	—	—	A	B
Molybdenum disul- fide	—	A	—	—
Monoethanolamine	B	D	A	A
N-Methyl-2-Pyrroli- done (NMP) ¹	A	C	A	—
Naphtha	D	A	B	A
Natural gas	A	A	A*	A
Nickel salts	A	A	A	B
Nitric acid (dil)	B	B	A	A
Nitric acid (med)	C	A	A	A
Nitric acid (conc)	D	A	A	A
Nitrobenzene	D	B	A	B
Nitrogen oxides	D	D	A	—
Nitrous acid	—	—	A	B
Oils, animal	B	A	A	A
Oils, mineral	B	A	A	A
Oils, vegetable	B	A	A	A
Oleic acid	D	B	A	B
Oxalic acid, cold	B	A	A	B
Oxygen, gas	B	B	A*	A
Palmitic acid, 100 % in ether	D	A	A	B

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Perchloric acid	D	A	A	C
Perchloroethylene	D	A	A	B
Phenol (carbolic acid)	D	A	A	B
Phosphoric acid, 50 %	C	A	A	A
Phthalic acid	B	B	A	B
Plating solutions	D	A	A	—
Polyglycol	A	A	—	—
Potassium carbonate	—	A	—	B
Potassium chlorate	B	A	A	B
Potassium hydroxide (med)	B	D	A	B
Potassium hydroxide (conc)	C	D	A	B
Potassium iodide	—	A	A	A
Propanol (propyl alcohol)	A	A	A	A
Pyridine	D	D	A	A
Silicone fluids	C	A	A	A
Silicone oils	C	A	A	A
Silver nitrate	A	A	A	B
Soap solutions	A	A	A	A
Sodium bicarbonate	A	A	A	B
Sodium bisulfate	—	A	A	D
Sodium bisulfite	A	A	A	B
Sodium borate	A	A	A	B
Sodium carbonate	A	A	A	A
Sodium chlorate	C	A	A	B
Sodium chloride	A	A	A	C
Sodium ferrocyanide	—	A	A	B
Sodium hydrosulfite	—	—	A	—

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Sodium hydroxide (dil)	A	A	A	A
Sodium hydroxide, 25 %	B	A	A	B
Sodium hydroxide (conc)	—	A	A	C
Sodium hypochlorite, <5%	B	A	A	A
Sodium hypochlorite, >5%	B	A	A	C
Sodium nitrate	D	A	A	B
Sodium silicate	A	A	A	B
Sodium sulfide	A	A	A	C
Sodium sulfite	A	A	A	A
Steam, up to 40 psi	A	B	A*	A
Stearic acid	B	A	A	B
Styrene	D	A	A	A
Sulfuric acid (dil)	D	A	A	D
Sulfuric acid (med)	D	A	A	D
Sulfuric acid (conc)	D	A	A	C
Sulfurous acid	D	B	A	B
Tannic acid	B	A	A	B
Tanning liquors	—	—	A	A
Tartaric acid	A	A	A	C
Tin salts	B	—	A	—
Toluene (toluol)	D	A	A	A
Trichloroacetic acid	D	C	A	D
Trichloroethylene	D	A	A	B
Trisodium phosphate	—	A	A	B
Turpentine	D	A	A	A
Urea	B	—	A	A

List of chemical resistance

Fluid Source: holland® APPLIED TECHNOLOGIES www.hollandapt.com	Silicone (peroxide/ platinum- cured)	Viton®	PTFE (Teflon)	Stainless steel
Uric acid	—	—	A	B
Water, fresh	B	A	A	A
Water, salt	A	A	A	B
Xylene	D	A	A	A
Zinc chloride	A	A	A	D

¹ see: <https://www.lyondellbasell.com/globalassets/documents/chemicals-technical-literature/elastomer-compatibility-2688>

3.5 Table 3

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Acetaldehyde (ethanal), 40%	A
Acetamide	A
Acetamide, 50%	A
Acetic acid butyl ester (butyl acetate)	A
Acetic acid (glacial acetic acid), concentrated	C
Acetic acid ethyl ester (ethyl acetate)	C
Acetic acid methyl ester (methyl acetate)	C
Acetic acid, 50%	C
Acetic acid, 10%	A
Acetic acid, 5%	A
Acetic ester (ethyl acetate, acetic acid ethyl ester)	C
Acetone (dimethyl ketone, propanone)	C
Acetophenone	A
Acetyl cellulose (cellulose acetate)	A
Acetylacetone	A
Acetylene gas	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Acetylene tetrachloride (tetrachloroethane)	A
Acrylic acid ethyl ester (ethyl acrylate)	A
Adipic acid	A
Adipic acid, aqueous	A
Aircraft engine fuels (JP)	A
Alcohol	A
Allyl alcohol	A
Allyl alcohol (2-propen-1-ol)	C
Alum (aluminum potassium sulfate), aqueous	A
Aluminum acetate, aqueous	A
Aluminum ammonium sulfate	C
Aluminum chloride	A
Aluminum fluoride	C
Aluminum hydroxide	A
Aluminum nitrate	C
Aluminum oxide	A
Aluminum phosphate, aqueous	A
Aluminum potassium sulfate (alum), aqueous	A
Aluminum potassium sulfate, diluted	C
Aluminum potassium sulfate, saturated	C
Aluminum sulfate	C
Aluminum sulfate, saturated	C
Amino acids	A
Aminoacetic acid	A
Aminobenzene (aniline)	C
Ammonia solution	A
Ammonia, 10%	A
Ammonia, aqueous (ammonia solution)	A
Ammonium acetate	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Ammonium carbonate (sal volatile)	A
Ammonium chloride (salmiac)	C
Ammonium chloride, aqueous	A
Ammonium chloride, solid	A
Ammonium fluoride (Fluorammon)	A
Ammonium hydroxide	A
Ammonium hydroxide, 5%	A
Ammonium hydroxide, 30%	C
Ammonium nitrate (fertilizer)	A
Ammonium nitrate, 10%	A
Ammonium nitrate, saturated	A
Ammonium phosphate (fertilizer)	A
Ammonium salts of mineral acids	A
Ammonium sulfate (fertilizer)	A
Ammonium sulfate, 10%	A
Ammonium sulfate, saturated	A
Ammonium thiocyanate	A
Amyl acetate (pentyl acetate)	A
Amyl alcohol (pentanol, pentyl alcohol)	A
Amyl chloride (pentyl chloride, 1-chloropentane)	A
Aniline (aminobenzene)	C
Animal fats	A
Anisole (methoxybenzene, methyl phenyl ether)	C
Anone (cyclohexanone)	A
Antichlor (sodium thiosulfate, fixing salt)	C
Antifreezes (glycol-based)	A
Argon gas	A
Barium carbonate	A
Barium chloride	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Barium hydroxide	A
Barium salts of mineral acids	A
Barium sulfate (barite, heavy spar)	A
Barium sulfide	A
Basic aluminium acetate	A
Beef tallow	A
Beer	A
Beeswax	A
Benzaldehyde (benzoic aldehyde)	A
Benzene	C
Benzoic acid (sodium salt)	A
Benzoic acid, aqueous	A
Benzoic acid, saturated	A
Benzoic aldehyde (benzaldehyde)	A
Benzoyl chloride	A
Benzyl alcohol (phenylcarbinol)	A
Benzyl carbinol (phenethyl alcohol)	A
Benzyl chloride (d-chlorotoluene)	A
Benzylamine	A
Bitumen	A
Bleaching powder	A
Blue vitriol (copper sulfate)	A
Borax (sodium tetraborate)	A
Boric acid	A
Boric acid, 10%	A
Boric acid, aqueous	A
Brake fluid, ATE blue (glycol based)	A
Brine	A
Bromobenzene (monobromobenzene)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Bromomethane	A
Bromotrifluoromethane	A
Butadiene	A
Butanal (butyraldehyde)	C
Butane	A
Butane, gaseous	A
Butane, liquid	A
Butanediol	A
Butanoic acid (butyric acid, propanecarboxylic acid)	C
Butanol (butyl alcohol)	A
Butanol-1 (butyl alcohol)	A
Butanone (methyl ethyl ketone, MEK)	A
Butene	A
Butter	A
Butylamine	C
Butyl acetate (acetic acid butyl ester)	A
Butyl alcohol (butanol)	A
Butyl glycol	A
Butyl glycol ether	A
Butyl methyl ether	A
Butyl oleate	A
Butylene glycol (2-butoxyethanol, butyl cellosolve)	A
Butylene, liquid	A
Butyraldehyde (butanal)	C
Butyric acid (butanoic acid, propanecarboxylic acid)	C
Calcium arsenate	A
Calcium bicarbonate	A
Calcium bisulfite	A
Calcium carbonate (chalk)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Calcium chloride	A
Calcium chloride, 10%	A
Calcium chloride, aqueous	A
Calcium hydroxide (milk of lime)	A
Calcium hydroxide (milk of lime), aqueous	A
Calcium hydroxide (milk of lime), concentrated	A
Calcium hypochlorite	A
Calcium hypochlorite, aqueous	C
Calcium hypochlorite, saturated	A
Calcium phosphate	A
Calcium phosphate, aqueous	A
Calcium sulfate (gypsum)	A
Caliche solution (Chile saltpeter)	A
Caraway	A
Carbolineum	A
Carbon dioxide	A
Carbon dioxide, wet	A
Carbon disulfide	A
Carbon monoxide, dry	A
Carbon tetrachloride (tetrachloromethane)	C
Carbon tetrachloride (tetrachloromethane, tetra)	C
Carbonic acid	A
Caustic potash (potassium hydroxide, potash lye)	A
Caustic soda (sodium hydroxide, soda lye)	A
Caustic soda (sodium hydroxide, soda lye), 50%	A
Cellosolve (ethyl glycol)	A
Cellulose acetate (acetyl cellulose)	A
Chalk (calcium carbonate)	A
Chlorine gas, dry	C

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Chlorobenzene (monochlorobenzene)	A
Chlorobromomethane (methylene chlorobromide)	C
Chlorodifluoromethane	A
Chloroethane (ethyl chloride)	A
Chloroethanol (chloroethyl alcohol, ethylene chlorohydrin)	C
Chlorofluoromethane	A
Chloromethane	A
Chloropentafluoroethane	A
Chlorothene (trichloroethane, methyl chloroform)	C
Chlorotrifluoromethane	A
Chromic acid, 10%	C
Citric acid	C
Citric acid, 10%	C
Citric acid, 50%	A
Citric acid, aqueous	C
Citric acid, saturated	A
Citrus fruits, citrus juices	A
Citrus juices	A
Citrus oils	A
Clophen (polychlorinated biphenyl)	A
Clophen A60, 50%	A
Coal gas	A
Coal tar oil	A
Cocoa	A
Cocoa butter	A
Cod liver oil	A
Colza oil	A
Common salt (sodium chloride)	A
Cooking fats, cooking oils	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Copper chloride, aqueous	C
Copper nitrate, aqueous	A
Copper sulfate (blue vitriol)	A
Copper sulfate (blue vitriol), 10%	A
Copper(II) nitrate, aqueous	A
Copper(II) nitrate, saturated	A
Cottonseed oil	A
Crude oil	A
Cumene (isopropylbenzene)	A
Cyclanon	A
Cycloaliphatic hydrocarbons	A
Cyclohexane (hexamethylene, hexahydrobenzene)	A
Cyclohexanol (hexalin, anol)	A
Cyclohexanone (anone)	A
Decahydronaphthalene (decane)	A
Decalin (decahydronaphthalene, decane)	A
Decane	A
Decanol	A
Dextrin	A
Dextrin, aqueous	A
Dextrorotatory d-tartaric acid	C
Diacetone (diacetone alcohol, pyranton)	A
Diacetone alcohol (pyranton, Dial, DA)	A
Dibenzyl ether	A
Dibromotetrafluoromethane	A
Dibutyl phthalate (palatinol A)	A
Dichlorobenzene	A
Dichlorodifluoromethane	A
Dichloroethane (ethylene chloride)	C

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Dichlorofluoromethane	A
Dichlorohexafluorocyclobutane	A
Dichloromethane (methylene chloride)	C
Dichlorotetrafluoroethane	A
Dicothyl phthalate	A
Diesel fuel	A
Diesel oil	A
Diethanolamine (DEA)	A
Diethyl ether (ether)	A
Diethyl ketone (3-pentanone)	A
Diethylamine	A
Diethylene glycol	A
Difluorochloroethane	A
Difluoroethane	A
Difluoromethane	A
Difluorotetrachloroethane	A
Dihydroxysuccinic acid (tartaric acid)	C
Diisobutyl ketone	A
Diisooctyl phthalate (DOP)	A
Diisopropyl ketone	A
Dimethyl ether	A
Dimethyl ketone (acetone, propanone)	C
Dimethyl phthalate (DMP)	A
Dimethyl sulfoxide (DMSO)	A
Dimethylbenzene (xylene)	A
Dimethylformamide (DMF)	A
Dinitrogen monoxide (laughing gas)	A
Dioctyl phthalate (octal, DOP)	A
Dioxane	C

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Dioxane, 1,4-	C
Diphenyl ether (diphenyl oxide)	C
Diphenyl oxide (diphenyl ether)	C
Diphyl	A
Dipropyl ketone	A
Dipropylene glycol	A
Dishwashing detergent	A
Disodium phosphate	A
Epichlorohydrin	A
Epsom salts (magnesium sulfate)	A
Ethanal (acetaldehyde)	A
Ethanol (ethyl alcohol, spirit), 100%	A
Ethanol (ethyl alcohol, spirit), 40%	A
Ethanol (ethyl alcohol, spirit), 50%	A
Ethanol (ethyl alcohol, spirit), 96%	A
Ethanol, 50%	A
Ethanol, 96%	A
Ether (diethyl ether)	A
Ethyl acetate (acetic acid ethyl ester, acetic ester)	C
Ethyl acrylate (acrylic acid ethyl ester)	A
Ethyl alcohol	A
Ethyl alcohol (ethanol, spirit)	A
Ethyl chloride (chloroethane)	A
Ethyl ether (ether, diethyl ether)	A
Ethyl glycol (cellosolve)	A
Ethyl methyl ketone	A
Ethylbenzene	A
Ethylene	A
Ethylene chloride (dichloroethane)	C

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Ethylene chlorohydrin (chloroethanol)	C
Ethylene glycol (glycol, 1,2-ethanediol)	A
Ethylene oxide	A
Ethylenediamine (1,2-diaminoethane)	A
Ethylenediaminetetraacetic acid (EDTA)	A
Ethylhexanol (isooctyl alcohol, 2-ethyl-1-hexanol)	A
Fats, cooking oils	A
Fats, mineral	A
Fats, vegetable and animal	A
Fatty acids	A
Fatty alcohol sulfonates	A
Ferric chloride (iron(III) chloride)	C
Ferrous chloride (iron(II) chloride)	C
Firewaters (ethyl alcohol, spirit)	A
Fixing baths	A
Fixing salt (antichlor, sodium thiosulfate)	C
Flaxseed oil (linseed oil)	A
Fluorammon (ammonium fluoride)	A
Fluorochloroalkanes (CFCs)	A
Formaldehyde (formalin, methanal)	A
Formaldehyde (formalin, methanal), 10%	A
Formaldehyde (formalin, methanal), 30%	A
Formaldehyde (formalin, methanal), 40%	A
Formamide	A
Formic acid methyl ester (methyl formate)	A
Formic acid, 3%	A
Freon (Frigen) 11	C
Freon (Frigen) 113	C
Fructose	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Fruit juices	A
Fruit pulp	A
Fruit wine	A
Fuel + 20% ethanol	A
Fuel + 20% methanol	A
Fuel oil, black coal- and brown coal-based	A
Fuel oil, mineral oil-based	A
Fuel, premium	A
Fuel, regular	A
Furfurol (furfuryl alcohol)	A
Gasoline	A
Gasoline, premium	A
Gasoline, unleaded	A
Gear oil, mineral	A
Gelatine	A
Gelatine, aqueous	A
Genantin (ethylene glycol)	A
Glacial acetic acid (acetic acid), concentrated	C
Glauber's salt (sodium sulfate)	A
Glucose (grape sugar)	A
Glue (animal glue)	A
Glycol (ethylene glycol)	A
Glysantin	A
Grape sugar (glucose)	A
Gum tragacanth	A
Glycerine (glycerol)	C
Gypsum (calcium sulfate)	A
HD Oil, motor oil free of aromatics	A
Helium gas	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Heptane	A
Hexachlorobenzene (HCB)	A
Hexachlorobutadiene (tripen)	A
Hexahydrobenzene (cyclohexane, hexamethylene)	A
Hexaldehyde	C
Hexalin (cyclohexanol)	A
Hexamethylene (cyclohexane, hexahydrobenzene)	A
Hexane	A
Hexanol (hexyl alcohol)	A
Hexyl alcohol (hexanol)	A
Honey	A
Hydraulic fluid, group H, HL, HLP (mineral oil-based)	A
Hydraulic fluid, group HSA (oil-in-water emulsion)	A
Hydraulic fluid, group HSB (water-in-oil emulsion)	A
Hydraulic fluid, group HSC (polyglycol-based)	A
Hydraulic fluid, group HSD a (R) (phosphoric ester)	A
Hydraulic fluid, group HSD b (S) (phosphoric ester)	A
Hydraulic fluid, group HSD c (T) (phosphoric ester)	A
Hydraulic oils (mineral oil-based)	A
Hydrazine	A
Hydrazine hydrate (diamide hydrate)	A
Hydrocyanic acid (hydrogen cyanide)	C
Hydrofluorosilicic acid (silicic acid)	A
Hydrogen peroxide, 3%	A
Hydrogen peroxide, 30%	A
Hydrogen peroxide, 90%	C
Hydrogen peroxide, diluted	A
Hydrogen sulfide	A
Hydrogen sulfide, saturated	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Hydrogen, gaseous	A
Ink	A
Iodine-potassium iodide solution	C
Iodine tincture	A
Iron(II) chloride, aqueous (ferrous chloride)	C
Iron(III) chloride (ferric chloride)	C
Iron nitrate, aqueous	A
Iron sulfate (iron vitriol), aqueous	A
Isoamyl alcohol	A
Isobutanol (isobutyl alcohol)	A
Isobutyl acetate	A
Isobutyl alcohol (isobutanol)	A
Isooctane	A
Isooctyl alcohol (ethylhexanol)	A
Isopropanol (2-propanol)	A
Isopropyl acetate	A
Isopropyl alcohol (isopropanol, persprit)	A
Isopropylacetone (methyl isobutyl ketone, MIBK)	A
Isopropylbenzene (cumene)	A
Kerosene	A
Lactic acid	A
Lactic acid, 3%	A
Lactic acid, 80%	C
Lactic acid, aqueous	C
Lactic acid, concentrated	C
Lactose	A
Lanolin (wool grease)	A
Latex	A
Laughing gas (nitrous oxide)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Lavender oil	A
Lead acetate (lead sugar)	A
Lead sugar (lead acetate)	A
Lead(II) acetate	A
Lead(II) acetate, aqueous	A
Lemon juice	A
Linseed oil (flaxseed oil)	A
Liqueurs	A
Liquid manure	A
Liquid soaps	A
Lithium bromide, aqueous	A
Lithium chloride, aqueous	A
Lubricating oil (mineral oil-based)	A
Lyes, diluted	A
Machine oil	A
Magnesium bromide	A
Magnesium carbonate	A
Magnesium chloride, aqueous	A
Magnesium chlorite	A
Magnesium hydroxide	A
Magnesium iodide	A
Magnesium nitrate	A
Magnesium sulfate (Epsom salt)	A
Margarine	A
MEK (butanone, methyl ethyl ketone)	A
Mercury	A
Mercury (II) chloride	C
Mercury chloride	C
Mercury chloride, aqueous	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Methanal (formaldehyde)	A
Methane	A
Methane (mine gas, natural gas)	A
Methanol (methyl alcohol, wood alcohol)	A
Methoxybenzene (anisole, methyl phenyl ether)	C
Methyl acetate (acetic acid methyl ester)	C
Methyl alcohol (methanol, wood alcohol)	A
Methyl butanol	A
Methyl butyl ether	C
Methyl cellosolve (methyl glycol acetate)	C
Methyl chloroform (trichloroethane, chloroethene)	C
Methyl formate (formic acid methyl ester)	A
Methyl glycol	A
Methyl glycol acetate (methyl cellosolve)	C
Methyl isobutyl ketone (isopropylacetone, MIBK)	A
Methylamine (monomethylamine)	A
Methylbenzene (toluene)	A
Methyl phenyl ether (anisole, methoxybenzene)	C
Methylene chloride (dichloromethane)	C
Methylene chlorobromide (chlorobromomethane)	C
Milk	A
Milk of lime (calcium hydroxide)	A
Mine gas (methane, natural gas)	A
Mineral oils	A
Mineral water	A
Molasses	A
Monobromobenzene	A
Monochloroacetic acid ester	A
Monochlorobenzene (chlorobenzene)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Mustard	A
Myristyl alcohol (myristic alcohol)	A
n-amyl acetate	A
Naphtha	A
Naphthalene (in alcohol)	A
Naphthalene (rock oil)	A
Natural gas (mine gas, methane)	A
Neon gas	A
Nickel nitrate	A
Nickel sulfate	A
Nickel sulfate, aqueous	A
Nickel sulfate, saturated	A
Nickel(II) chloride	A
Nitrobenzene (oil of mirbane)	C
Nitrogen	A
Nicotine, aqueous	A
Nicotinic acid	A
Nitrotoluene	C
Octafluorocyclobutane	A
Octal (dioctyl phthalate, DOP)	A
Octane	A
Octyl alcohol	A
Oil of mirbane (nitrobenzene)	C
Oleic acid	A
Olive oil	A
Orange peel oil	A
Oxalic acid (ethanedioic acid)	A
Oxygen	A
P3 solution, aqueous	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Palatinol A (dibutyl phthalate)	A
Palmitic acid	C
Paraffin (alkane)	A
Paraffin oil	A
Paraffin wax, molten	A
Peanut oil	A
Pectin	A
Pectin, aqueous	A
Pentane	A
Pentanol (pentyl alcohol, amyl alcohol)	A
Pentyl alcohol (pentanol, amyl alcohol)	A
Pentyl chloride (amyl chloride, 1-chloropentane)	A
Perchloric acid	C
Perchloric acid (hyperchloric acid)	C
Perchloroethylene (tetrachloroethylene)	C
Perfluoropropane	A
Petroleum	A
Petroleum ether	A
Phenylcarbinol (benzyl alcohol)	A
Phenylethyl alcohol (benzyl carbinol)	A
Phenylhydrazine	A
Phosgene, liquid	C
Phosphine	A
Phosphoric acid, 1-5%	C
Phosphoric acid, 10%	C
Phosphoric acid, 85%	A
Phosphorus trichloride	C
Photo developers	A
Photographic developers	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Photographic emulsions	A
Phthalic acid (benzene-o-dicarboxylic acid)	A
Picric acid, aqueous	C
Polyethylene glycol	A
Polysolvan O	A
Potassium acetate	A
Potassium bichromate, 10%	A
Potassium borate, 10%	A
Potassium bromide	A
Potassium carbonate	A
Potassium carbonate (potash)	A
Potassium carbonate (potash), aqueous	A
Potassium carbonate, aqueous	A
Potassium chlorate, aqueous	A
Potassium chlorate, saturated	A
Potassium chloride (sylvine)	A
Potassium cyanide	A
Potassium cyanide, aqueous	C
Potassium cyanide, saturated	C
Potassium dichromate, aqueous	C
Potassium dichromate, saturated	C
Potassium ferricyanide	A
Potassium ferrocyanide	A
Potassium fluoride	A
Potassium hydrogen carbonate	A
Potassium hydrogen sulfate	A
Potassium hydrogen sulfate, aqueous	A
Potassium hydroxide (caustic potash, potash lye)	A
Potassium hydroxide (caustic potash, potash lye), 1%	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Potassium hydroxide (caustic potash, potash lye), 10%	A
Potassium hydroxide (caustic potash, potash lye), 50%	C
Potassium iodide, aqueous	A
Potassium iodide, saturated	A
Potassium nitrate (saltpeter)	A
Potassium nitrate (saltpeter), 50%	A
Potassium nitrate (saltpeter), aqueous	A
Potassium perchlorate, aqueous	A
Potassium perchlorate, saturated	A
Potassium permanganate	C
Potassium permanganate, 1%	A
Potassium permanganate, aqueous	C
Potassium persulfate	C
Potassium sulfate	C
Potassium sulfite	A
Potassium thiosulfate	A
Propane, gaseous	A
Propanecarboxylic acid (butyric acid, butanoic acid)	C
Propane, liquid	A
Propanediol (propylene glycol)	A
Propanol (propyl alcohol)	A
Propanone (acetone, dimethyl ketone)	C
Propene	A
Propyl acetate	A
Propyl alcohol (propanol)	A
Propylamine	C
Propylene	A
Propylene glycol (propanediol)	A
Prussic acid (hydrocyanic acid)	C

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Pyranton (diacetone alcohol, Dial, DA)	A
Pyridine	C
Pyrogallol (pyrogalllic acid)	C
Rubber dispersion	A
Sal ammoniac (ammonium chloride)	C
Sal volatile (ammonium carbonate)	A
Salt water (sea water)	A
Saltpeter (potassium nitrate)	A
Sea water (salt water)	A
Silicic acid (hydrofluorosilicic acid)	A
Silicone grease	A
Silicone oil	A
Silver acetate	C
Silver nitrate (lunar caustic)	C
Silver nitrate (lunar caustic), aqueous	A
Soap water (soap hydrous solution)	A
Soda (sodium carbonate)	A
Soda (sodium carbonate), aqueous	A
Sodium acetate	A
Sodium benzoate	A
Sodium benzoate, aqueous	A
Sodium bicarbonate, aqueous	A
Sodium bisulfate	A
Sodium bisulfate, aqueous	A
Sodium borate, aqueous	A
Sodium borate, saturated	A
Sodium bromate	A
Sodium bromide	A
Sodium carbonate (soda)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Sodium carbonate (soda), aqueous	A
Sodium chlorate	A
Sodium chlorate, aqueous	A
Sodium chloride (common salt)	A
Sodium cyanide	A
Sodium cyanide, aqueous	C
Sodium cyanide, saturated	C
Sodium dichromate	A
Sodium dodecylbenzenesulfonate	A
Sodium fluoride	A
Sodium hexametaphosphate	A
Sodium hydroxide (caustic soda, soda lye)	A
Sodium hydroxide (caustic soda, soda lye), 1%	A
Sodium hydroxide (caustic soda, soda lye), 30%	A
Sodium hydroxide (caustic soda, soda lye), 45%	A
Sodium hydroxide (caustic soda, soda lye), 50%	A
Sodium iodide	A
Sodium nitrate (Chile saltpeter)	A
Sodium nitrate (Chile saltpeter), aqueous	A
Sodium nitrate (Chile saltpeter), saturated	A
Sodium nitrite, aqueous	A
Sodium nitrite, saturated	A
Sodium oxalate	A
Sodium perborate	A
Sodium perborate, aqueous	A
Sodium perborate, saturated	A
Sodium perchlorate	A
Sodium phosphate	A
Sodium phosphate, aqueous	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Sodium phosphate, saturated	A
Sodium silicate, aqueous	A
Sodium stearate	A
Sodium sulfate (Glauber's salt)	A
Sodium sulfate (Glauber's salt), aqueous	A
Sodium sulfate (Glauber's salt), saturated	A
Sodium sulfide	A
Sodium sulfide, aqueous	A
Sodium sulfide, saturated	A
Sodium sulfite	A
Sodium tetraborate (borax)	A
Sodium thiosulfate (antichlor, fixing salt)	C
Soft soap	A
Soybean oil	A
Spermaceti	A
Spirit (ethyl alcohol, ethanol)	A
Spruce oil (spruce needle oil)	A
Starch	A
Starch sirup	A
Starch, aqueous	A
Stearic acid (zinc salt)	A
Styrene (styrol)	A
Styrene (vinyl benzene, phenylethene)	A
Succinic acid	C
Suds	A
Sugar beet juice	A
Sugar sirup	A
Sulfur	A
Sulfur hexafluoride	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Tallow	A
Tannic acid (tannin)	C
Tannin (tannic acid)	C
Tar	A
Tartaric acid (dihydroxysuccinic acid)	C
Tetrachloroethane (acetylene tetrachloride)	A
Tetrachloroethylene (perchloroethylene)	C
Tetrachloromethane (carbon tetrachloride)	C
Tetrafluoromethane	A
Tetrahydrofuran (diethylene oxide, tetramethylene oxide)	C
Tetrahydronaphthalene (tetralin)	A
Tetralin (tetrahydronaphthalene)	A
Tetramethylene oxide (tetrahydrofuran)	C
Thiofuran (thiol, thiophene)	C
Thionyl chloride	A
Thiophene (thiofuran, thiol)	C
Toluene (methylbenzene)	A
Transformer oil	A
Transformer oil (mineral oil-based)	C
Trichloroethane (chloroethene, methyl chloroform)	C
Trichlorofluoromethane	A
Trichlorotrifluoroethane	A
Tricresyl phosphate	A
Triethanolamine (triethylolamide)	A
Triethylamine (TEA)	A
Triethylene glycol	A
Trifluorotrichloroethane	A
Triglycol	A
Triglycol acetate	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Trilone B, 10%	A
Trimethylolpropane	A
Tripfen (hexachlorobutadiene)	A
Tripropylene glycol (TPG)	A
Trisodium phosphate	A
Turbine oil, mineral	A
Turpentine	A
Two-stroke (engine) oil	A
Undecyl alcohol	A
Urea	A
Urea, aqueous	A
Uric acid	A
Urine	A
Vaseline	A
Vaseline oil	A
Vegetable oil	A
Vinegar	A
Water, distilled	A
Water, with carbon dioxide	A
Water/Spring water, cold	A
Wax molten	A
Whisky	A
White spirit	A
Wine, brandy	A
Wood alcohol (methanol, methyl alcohol)	A
Wood oil	A
Wool grease (lanolin)	A
Xenon	A
Xylene (dimethylbenzene)	A

List of chemical resistance

Fluid Source: Reichelt Chemietechnik https://www.rct-online.de	POM
Yeast	A
Zinc carbonate	A
Zinc chloride	A
Zinc chloride, 10%	A
Zinc chloride, aqueous	C
Zinc oxide	A
Zinc phosphate	A
Zinc salt (stearic acid)	A
Zinc stearate	A
Zinc sulfate	A
Zinc sulfate, 10%	A

The CHEMICAL PROPERTIES of N-BK7 glass see: https://www.schott.com/d/advanced_optics/c36214d9-13c4-468c-bf40-8d438b89f532/1.30/schott-optical-glass-pocket-catalog-jan-2020-row.pdf



